

m|devices®

AN140002NS - AN140003NS Anaesthetic Masks

Scented

INSTRUCTION FOR USE

DESCRIPTION

m|devices scented Anaesthetic Masks come in two color-coded sizes for infants and children, ensuring easy identification.

The soft facial cushion creates a secure seal for optimal delivery of oxygen or anesthetic gas.

The one-way valve supports efficient and easy inflation and deflation.

FOR USE BY A QUALIFIED CLINICIAN. THE BELOW IS ONLY A SUGGESTION AND FACILITY PROTOCOL MUST BE FOLLOWED FOR ALL CLINICAL PROCEDURES WHERE THIS PRODUCT IS USED.

CAUTION

- Single use only.
- DO NOT re-sterilise.
- DO NOT store at extreme temperatures and humidity, avoid direct sunlight. Handle with care.
- NON STERILE, DO NOT use if the package or product has been damaged or contaminated.
- EU Notice: any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and /or patient is established.

STEPS

1. Inspect cushion for inflation pressure.
2. If required, use a syringe to inflate or deflate air cushion through one-way valve.
3. Connect mask to breathing circuit or resuscitation device.
4. Ensure connections are firm.
5. Cover the patients face with mask and apply gentle but firm pressure to create an effective mask seal.

PRECAUTION


- Test the Anaesthetic Mask on the patient's face first to ensure cushion is adequately inflated to obtain an adequate seal prior to connecting to breathing circuit or resuscitation device.

m|devices

27 Llewellyn Avenue
Villawood NSW 2163
Australia
T +61(0)2 8718 2800
F +61(0)2 8718 2801

 AIRWAY
MANAGEMENT Anaesthesia Masks

 **Guangzhou Shengwei
Medical Devices, CO., LTD**
No. 126, Guangshao Road
Aotou Town, Conghua District
Guangzhou 510945
Guangdong Province, P.R. China

 **MedPath GmbH**
Mies-van-der-Rohe-Strasse 8
80807 Munich
Germany



Rev.05