

m|devices®

IV052416

Standard Bore Extension Set

Needleless Access Site

Needleless to Male Luer Lock with Rotating Collar

15cm



## INSTRUCTION FOR USE

### DESCRIPTION

The m|devices IV Standard Bore Extension Set is a single-lumen extension set with needleless access site. It is 15cm in length with a male Luer lock and rotating collar to ensure a secure connection for your IV administration.

**FOR USE BY A QUALIFIED CLINICIAN. THE BELOW IS ONLY A SUGGESTION AND FACILITY PROTOCOL MUST BE FOLLOWED FOR ALL CLINICAL PROCEDURES WHERE THIS PRODUCT IS USED.**

### CAUTION

- Single use only.
- DO NOT re-sterilise.
- DO NOT store at extreme temperatures and humidity, avoid direct sunlight. Handle with care.
- STERILE (EO), DO NOT use if the package or product has been damaged or contaminated.
- EU Notice: any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and /or patient is established.

### STEPS

Ensure aseptic technique is followed at all times.

1. Check that all connection components are securely fastened.
2. Prime the line by inserting a pre-filled syringe into the Needleless Access Site and depressing the plunger to flush fluid through the set, ensuring all air bubbles are removed.
3. Alternatively, connect to an existing IV administration set by closing the roller regulator, connecting the distal end of the administration set, and opening the roller regulator to allow fluid to flow until it drips out of the distal end of the set.
4. Remove the protective cap from the male Luer connector.
5. Connect and secure it to the access device or peripheral catheter by inserting it and engaging the rotating collar.


### PRECAUTIONS


- The security of all connections should be checked once established and monitored during use to prevent disconnection.

m|devices

27 Llewellyn Avenue  
Villawood NSW 2163  
Australia  
T +61(0)2 8718 2800  
F +61(0)2 8718 2801



 **BQ Plus Medical Co., Ltd.**  
No. 18, Che Ye Road  
Che Dun Town, Songjiang  
201611 Shanghai, China

 **Prolinx GmbH**  
Brehmstr. 56  
40239 Duesseldorf  
Germany



Rev New