




CONTINENCE  
CARE

Urology &  
Continence  
Products

m|devices®







We believe in effortless performance, so you can focus on caregiving.

One size never fits all. For over a decade we've been designing and manufacturing devices that simplify the roles of clinicians and carers – not just cut to a standard specification but tailored to the unique needs of Australians.

Our focus is on making performance effortless, leaving you free to focus on what matters most – caring for patients.



Leading by Design is a philosophy that reflects our commitment to research and collaboration. A crucial part of our success has been our effort to embed medical professionals in the design process. Their insight helps inform solutions that make the giving of care simple and rewarding.

The problems you face today are the solutions we're working towards for tomorrow.





## CONTINENCE CARE

### Accessible, reliable and simple to use

The treatment of acute and chronic urinary disorders affects many groups in our community and is often an invasive and confronting process for the patient.

As care givers and clinicians, you need the reliability and performance that comes from solutions carefully tailored to your needs, enabling you to perform at your best and assure comfort and confidence for those in your care.

Our range of catheters and accessories are designed to simplify the care giving process, putting you in control of a reassuring, patient-centred experience that minimises discomfort.

### Contents

#### Catheters

##### Indwelling Catheters

- 07 | 2-Way Foleys
- 08 | 2-Way Foleys (Open-Ended)
- 09 | 2-Way Foleys (Tiemann Tip)
- 10 | 3-Way Foleys

##### Intermittent Catheters

- 12 | Nelaton Catheters
- 13 | Nelaton Catheters (Hydrophilic Coated)
- 14 | Nelaton Catheters (Hydrophilic with Water Sachet)
- 15 | Nelaton Catheters (Hydrophilic with Water Sachet - Coude Tip)
- 16 | Intermittent Closed Catheterisation System

#### Drainage Bags

##### Urine Bags

- 19 | Urine Bags (Outlet)
- 19 | Urine Bags (No Outlet)

##### Leg Bags

- 20 | Leg Bags (Lever Tap)
- 21 | Leg Bags (T-Tap)

#### Accessories

- 23 | Universal Securement Device
- 24 | Urine Bag Hanger
- 25 | Catheter Tip Syringe






## Catheters

Our catheters are designed to support your full range of acute and chronic care needs.


Available in paediatric, female and male sizes, each colour coded product has been developed to simplify product selection, maximise patient comfort and support an intuitive, confident insertion.

### 2-Way Foley Catheters


Approved for both urethral and suprapubic\* catheterisation. They feature a symmetrical balloon design to support reliable positioning and the unisex range comes with high contrast incremental markings to aid catheter placement.



- ⊗ Intended for Urethral and SPC use
- ⊗ 100% Silicone
- ⊗ High contrast incremental markings from 6cm to 30cm
- ⊗ Long indwelling time of up to 8 weeks\*\*
- ⊗ Smooth rounded tip to aid gentle insertion
- ⊗ Symmetrical balloon ensures reliable positioning
- ⊗ Three lengths available: Paediatric 33cm, Female 23cm and Unisex 45cm
- ⊗ Colour coded for easy size identification



Round Tip



Incremental markings to support placement

#### 2-Way Foley Catheters

RE-ORDER	DESCRIPTION	GAUGE	COLOUR	STERILE CARTON QTY
UR010006	2-Way Foley Catheter, Standard Tip, 33cm (Paediatric) with 1.5mL Balloon	6Fr	Pink	10
UR010000	2-Way Foley Catheter, Standard Tip, 33cm (Paediatric) with 3mL Balloon	8Fr	Light Blue	10
UR010001	2-Way Foley Catheter, Standard Tip, 33cm (Paediatric) with 3mL Balloon	10Fr	Black	10
UR010012	2-Way Foley Catheter, Standard Tip, 23cm (Female) with 10mL Balloon	12Fr	White	10
UR010014	2-Way Foley Catheter, Standard Tip, 23cm (Female) with 10mL Balloon	14Fr	Green	10
UR010016	2-Way Foley Catheter, Standard Tip, 23cm (Female) with 10mL Balloon	16Fr	Orange	10
UR011000	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	12Fr	White	10
UR011001	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	14Fr	Green	10
UR011002	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	16Fr	Orange	10
UR011003	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	18Fr	Red	10
UR011004	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	20Fr	Yellow	10
UR011005	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	22Fr	Purple	10
UR011006	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	24Fr	Blue	10
UR011007	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	26Fr	Pink	10
UR011016	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	16Fr	Orange	10
UR011017	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	18Fr	Red	10
UR011018	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	20Fr	Yellow	10
UR011019	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	22Fr	Purple	10
UR011020	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	24Fr	Blue	10
UR011021	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	26Fr	Pink	10

\* SPC use excludes paediatric range.


\*\* As per laboratory-based testing results and validation against requirements ISO 20696.



## 2-Way Foley Catheters (Open-Ended)

With its open-ended design and bilateral eyelets this Open Ended Foley Catheter is designed to deliver the high flow rates for urine drainage necessary to support patients with a strong tendency for catheter blockages.


It features clear measurements to support catheter placement.



- ⊗ Incremental markings to support catheter placement
- ⊗ 100% Silicone
- ⊗ Short tip to reduce bladder irritation
- ⊗ Open ended with two bilateral eyelets to improve flow
- ⊗ Suitable for both males and females
- ⊗ Symmetrical balloon ensures reliable positioning
- ⊗ Intended for urethral and SPC use
- ⊗ Long indwelling time of up to 8 weeks\*
- ⊗ Colour coded for easy size identification
- ⊗ Radiopaque



Open Tip to improve flow



Incremental markings to support placement

### 2-Way Foley Catheters (Open-Ended) STERILE


RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
UR074014	Open-Ended 2-Way Foley Catheter, 40cm with 10mL Balloon	14Fr	Green	10
UR074016	Open-Ended 2-Way Foley Catheter, 40cm with 10mL Balloon	16Fr	Orange	10
UR074018	Open-Ended 2-Way Foley Catheter, 40cm with 10mL Balloon	18Fr	Red	10
UR074020	Open-Ended 2-Way Foley Catheter, 40cm with 10mL Balloon	20Fr	Yellow	10
UR074022	Open-Ended 2-Way Foley Catheter, 40cm with 10mL Balloon	22Fr	Purple	10
UR074024	Open-Ended 2-Way Foley Catheter, 40cm with 10mL Balloon	24Fr	Blue	10

\* As per laboratory-based testing results and validation against requirements ISO 20696.

## 2-Way Foley Catheters (Tiemann Tip)

Equipped with all the features of traditional Foley catheters the Tiemann Tip version is designed to support the more demanding insertion.

- ⊗ 100% Silicone
- ⊗ Tiemann tip supports more demanding insertions
- ⊗ Colour-coded to support easy size identification
- ⊗ Symmetrical balloon ensures reliable positioning
- ⊗ Radiopaque



Tiemann Tip

### 2-Way Foley Catheters (Tiemann Tip) STERILE


RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
UR012000*	2-Way Foley Catheter, Tiemann Tip, 41cm with 10mL Balloon	12Fr	White	10
UR012001*	2-Way Foley Catheter, Tiemann Tip, 41cm with 10mL Balloon	14Fr	Green	10
UR012002*	2-Way Foley Catheter, Tiemann Tip, 41cm with 10mL Balloon	16Fr	Orange	10
UR012003*	2-Way Foley Catheter, Tiemann Tip, 41cm with 10mL Balloon	18Fr	Red	10
UR012004*	2-Way Foley Catheter, Tiemann Tip, 41cm with 10mL Balloon	20Fr	Yellow	10
UR012005*	2-Way Foley Catheter, Tiemann Tip, 41cm with 10mL Balloon	22Fr	Purple	10
UR012006*	2-Way Foley Catheter, Tiemann Tip, 41cm with 10mL Balloon	24Fr	Blue	10
UR012007*	2-Way Foley Catheter, Tiemann Tip, 41cm with 10mL Balloon	26Fr	Pink	10



## 3-Way Foley Catheters

The 3-Way Foley is made from 100% medical grade translucent and flexible silicone. The symmetrical balloon design supports reliable positioning in the bladder and minimises the potential for ridging.

- ⊗ 100% Silicone
- ⊗ Extended 43cm length
- ⊗ Long indwelling time of up to 29 days
- ⊗ Colour-coded to support easy size identification
- ⊗ Symmetrical balloon ensures reliable positioning
- ⊗ Radiopaque



3-Way Foley Catheters					STERILE
RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY	
UR013000	3-Way Foley Catheter, Standard Tip, 40cm with 30mL Balloon	16Fr	Orange	10	
UR013001	3-Way Foley Catheter, Standard Tip, 43cm with 30mL Balloon	18Fr	Red	10	
UR013002	3-Way Foley Catheter, Standard Tip, 43cm with 50mL Balloon	20Fr	Yellow	10	
UR013003	3-Way Foley Catheter, Standard Tip, 43cm with 50mL Balloon	22Fr	Purple	10	
UR013004	3-Way Foley Catheter, Standard Tip, 43cm with 60mL Balloon	24Fr	Blue	10	

## Where confidence meets care.

Our focus is on making performance effortless, leaving you free to focus on what matters most – caring for patients.






## Nelaton Catheters

Our Nelaton Catheters are intended for intermittent bladder catheterisation. They are made from medical grade thermo-sensitive plastic and offer a perfectly finished distal end for smooth insertion. The blue coloured holding sleeve assists to maintain aseptic technique during insertion.

- ⊗ Available in male, female and paediatric lengths
- ⊗ Smooth distal end aids insertion
- ⊗ Colour-coded to support easy size identification
- ⊗ Holding sleeve to maintain aseptic technique during insertion




Nelaton Catheters				STERILE
RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
UR040000*	Standard Nelaton Catheter, 30cm (Paediatric)	6Fr	Light Green	50
UR040001*	Standard Nelaton Catheter, 30cm (Paediatric)	8Fr	Light Blue	50
UR040002	Standard Nelaton Catheter, 30cm (Paediatric)	10Fr	Black	50
UR041000	Standard Nelaton Catheter, 20cm (Female)	8Fr	Light Blue	50
UR041001	Standard Nelaton Catheter, 20cm (Female)	10Fr	Black	50
UR041002	Standard Nelaton Catheter, 20cm (Female)	12Fr	White	50
UR041003	Standard Nelaton Catheter, 20cm (Female)	14Fr	Green	50
UR041004	Standard Nelaton Catheter, 20cm (Female)	16Fr	Orange	50
UR041005	Standard Nelaton Catheter, 20cm (Female)	18Fr	Red	50
UR041006	Standard Nelaton Catheter, 20cm (Female)	20Fr	Yellow	50
UR042000	Standard Nelaton Catheter, 40cm (Male)	8Fr	Light Blue	50
UR042001	Standard Nelaton Catheter, 40cm (Male)	10Fr	Black	50
UR042002	Standard Nelaton Catheter, 40cm (Male)	12Fr	White	50
UR042003	Standard Nelaton Catheter, 40cm (Male)	14Fr	Green	50
UR042004	Standard Nelaton Catheter, 40cm (Male)	16Fr	Orange	50
UR042005	Standard Nelaton Catheter, 40cm (Male)	18Fr	Red	50
UR042006*	Standard Nelaton Catheter, 40cm (Male)	20Fr	Yellow	50

\* This product is currently being sold without the blue sleeve but will transition to blue sleeve upon depletion of existing stock.

## Nelaton Catheters (Hydrophilic Coated)

Our Hydrophilic catheters are made from non-toxic DEHP free PVC. Their smooth water soluble polymer coating combined with a heat polished distal tip supports a smooth insertion without the need of lubricants.


- ⊗ Available in male, female and paediatric lengths
- ⊗ Heat polished distal end aids insertion
- ⊗ Hydrophilic coating supports lubricant free insertion
- ⊗ Colour coded to support easy size identification
- ⊗ DEHP free



Nelaton Catheters (Hydrophilic Coated)				STERILE
RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
UR043001	Hydrophilic Coated Nelaton Catheter, 30cm (Paediatric)	8Fr	Light Blue	25
UR043002	Hydrophilic Coated Nelaton Catheter, 30cm (Paediatric)	10Fr	Black	25
UR044000	Hydrophilic Coated Nelaton Catheter, 18cm (Female)	8Fr	Light Blue	25
UR044001	Hydrophilic Coated Nelaton Catheter, 18cm (Female)	10Fr	Black	25
UR044002	Hydrophilic Coated Nelaton Catheter, 18cm (Female)	12Fr	White	25
UR044003	Hydrophilic Coated Nelaton Catheter, 18cm (Female)	14Fr	Green	25
UR044004	Hydrophilic Coated Nelaton Catheter, 18cm (Female)	16Fr	Orange	25
UR044005	Hydrophilic Coated Nelaton Catheter, 18cm (Female)	18Fr	Red	25
UR045000	Hydrophilic Coated Nelaton Catheter, 40cm (Male)	8Fr	Light Blue	25
UR045001	Hydrophilic Coated Nelaton Catheter, 40cm (Male)	10Fr	Black	25
UR045002	Hydrophilic Coated Nelaton Catheter, 40cm (Male)	12Fr	White	25
UR045003	Hydrophilic Coated Nelaton Catheter, 40cm (Male)	14Fr	Green	25
UR045004	Hydrophilic Coated Nelaton Catheter, 40cm (Male)	16Fr	Orange	25
UR045005	Hydrophilic Coated Nelaton Catheter, 40cm (Male)	18Fr	Red	25

## Nelaton Catheters (Hydrophilic with Water Sachet)

Our Hydrophilic Nelaton Catheter with Water Sachet is an “all-in-one” product which is convenient and easy to use. It is available in various lengths and comes with a holding sleeve.



- ⊗ No-touch sleeve
- ⊗ Easy peel packaging with wall attachment
- ⊗ Hydrophilic coating supports lubricant free insertion
- ⊗ Heat polished distal end aids insertion
- ⊗ Available in male, female and paediatric lengths


### Nelaton Catheters (Hydrophilic with Water Sachet)

RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY	STERILE
UR046000	Hydrophilic Nelaton Catheter with Water Sachet, 30cm (Paediatric)	6Fr	Light Green	25	
UR046001	Hydrophilic Nelaton Catheter with Water Sachet, 30cm (Paediatric)	8Fr	Light Blue	25	
UR046002	Hydrophilic Nelaton Catheter with Water Sachet, 30cm (Paediatric)	10Fr	Black	25	
UR047000	Hydrophilic Nelaton Catheter with Water Sachet, 18cm (Female)	8Fr	Light Blue	25	
UR047001	Hydrophilic Nelaton Catheter with Water Sachet, 18cm (Female)	10Fr	Black	25	
UR047002	Hydrophilic Nelaton Catheter with Water Sachet, 18cm (Female)	12Fr	White	25	
UR047003	Hydrophilic Nelaton Catheter with Water Sachet, 18cm (Female)	14Fr	Green	25	
UR047004	Hydrophilic Nelaton Catheter with Water Sachet, 18cm (Female)	16Fr	Orange	25	
UR047005	Hydrophilic Nelaton Catheter with Water Sachet, 18cm (Female)	18Fr	Red	25	
UR048000	Hydrophilic Nelaton Catheter with Water Sachet, 40cm (Male)	8Fr	Light Blue	25	
UR048001	Hydrophilic Nelaton Catheter with Water Sachet, 40cm (Male)	10Fr	Black	25	
UR048002	Hydrophilic Nelaton Catheter with Water Sachet, 40cm (Male)	12Fr	White	25	
UR048003	Hydrophilic Nelaton Catheter with Water Sachet, 40cm (Male)	14Fr	Green	25	
UR048004	Hydrophilic Nelaton Catheter with Water Sachet, 40cm (Male)	16Fr	Orange	25	
UR048005	Hydrophilic Nelaton Catheter with Water Sachet, 40cm (Male)	18Fr	Red	25	

## Nelaton Catheters (Hydrophilic with Water Sachet - Coude Tip)

Our Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip is an “all-in-one” product which is convenient and easy to use.

The angled tip makes it easy to navigate the urethra and supports a difficult insertion.



Coude Tip

- ⊗ Supports a difficult insertion
- ⊗ No touch sleeve
- ⊗ Easy peel packaging with wall attachment
- ⊗ Hydrophilic coating supports lubricant free insertion
- ⊗ Heat polished distal end aids insertion

### Nelaton Catheters (Hydrophilic with Water Sachet - Coude Tip)


RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY	STERILE
UR049001	Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip, 40cm (Male)	10Fr	Black	25	
UR049002	Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip, 40cm (Male)	12Fr	White	25	
UR049003	Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip, 40cm (Male)	14Fr	Green	25	
UR049004	Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip, 40cm (Male)	16Fr	Orange	25	
UR049005	Hydrophilic Nelaton Catheter with Water Sachet - Coude Tip, 40cm (Male)	18Fr	Red	25	



## Intermittent Closed Catheterisation System

Our Closed Catheterisation System is sterile and pre-lubricated with a hydrophilic gel, so no additional preparation is needed. The system works like a traditional Nelaton Catheter and provides a simple catheterisation that reduces the risk of urinary tract infections.

- ⊙ Integrated 1500mL bag
- ⊙ Sterile closed system reduces the chance of infection
- ⊙ Pre-lubricated for simple insertion
- ⊙ DEHP, BPA and latex free



Intermittent Closed Catheterisation System			STERILE
RE-ORDER	DESCRIPTION	GAUGE	CARTON QTY
UR050008	Standard Intermittent Catheter with Gel and 1500mL Bag	8Fr	25
UR050010	Standard Intermittent Catheter with Gel and 1500mL Bag	10Fr	25
UR050012	Standard Intermittent Catheter with Gel and 1500mL Bag	12Fr	25
UR050014	Standard Intermittent Catheter with Gel and 1500mL Bag	14Fr	25
UR050016	Standard Intermittent Catheter with Gel and 1500mL Bag	16Fr	25



Designed for medical professionals, by medical professionals.

We know from experience that by observing, listening and collaborating with the right people we can find performance breakthroughs with the potential to transform how you work.





## Urine Bags

Our Urine Bags are available with and without a non-return valve and variations in tube length ensure the bag is optimised for different settings.

- ⊗ Clear volume graduation to aid simple volume measurement
- ⊗ Varying tube lengths to support different patient needs
- ⊗ Made from latex-free, durable material to deliver confidence



## Drainage Bags

Our drainage bags are designed to be unobtrusive and minimise the risk of unwanted spills providing confidence and security for your patients.

We provide a range of different urine bags to suit lifestyle needs and the care context. We also offer a double and triple chamber integrated meter to simplify the monitoring of urine flow.

Urine Bags				STERILE
RE-ORDER	DESCRIPTION	TUBING LENGTH	OUTLET	CARTON QTY
UR062000	Urine Bag – 2000mL T-TAP NRV 90cm Tubing	90cm	T-Tap	250
UR062001	Urine Bag – 2000mL T-TAP 90cm Tubing	90cm	T-Tap	250
UR062003	Urine Bag – 2000mL T-TAP NRV 120cm Clear Tubing, Sample Port	120cm	T-Tap	100
UR062007	Urine Bag – 2000mL T-TAP NRV 110cm Tubing	110cm	T-Tap	250
UR062008	Urine Bag – 2000mL T-TAP 110cm Tubing	110cm	T-Tap	250
UR062500	Urine Bag – 2500mL NRV, Integrated Hanger 120cm Tubing	120cm	T-Tap	40
UR064000	Urine Bag – 4000mL T-TAP NRV Integrated with Rope 150cm Tubing	150cm	T-Tap	40
RE-ORDER	DESCRIPTION	TUBING LENGTH	OUTLET	CARTON QTY
UR062004	Urine Bag – 2000mL NRV Hanger 120cm Tubing (no outlet)	120cm	No outlet	50
UR062005	Urine Bag – 2000mL NRV 120cm Tubing (no outlet)	120cm	No outlet	100



## Leg Bags (Lever Tap)

Our Leg Bags are designed to be unobtrusive offering discreet storage under clothing. They have a Lever Tap design to support users with reduced dexterity. To improve comfort, they feature a soft non-woven backing, silicone lined strap and fabric fastener.

- Discreet to wear
- Soft non-woven backing is comfortable next to the skin
- Silicone lined straps and fabric fasteners
- Multiple lengths: 10cm, 30cm and 50cm
- Various volume capacities; 350mL, 500mL and 750mL
- Bonded step connector
- Support users who have reduced manual dexterity as the Lever Tap can easily be opened and closed



Leg Bags (Lever Tap)			STERILE
RE-ORDER	DESCRIPTION	CAPACITY	CARTON QTY
UR335010*	LEG BAG - Lever Tap NRV 10cm Short Tube with Silicone lined Straps	350mL	80
UR335030*	LEG BAG - Lever Tap NRV 30cm Long Tube with Silicone lined Straps	350mL	80
UR335050*	LEG BAG - Lever Tap NRV 50cm Long Tube with Silicone lined Straps	350mL	80
UR350010	LEG BAG - Lever Tap NRV 10cm Short Tube with Silicone lined Straps	500mL	80
UR350030*	LEG BAG - Lever Tap NRV 30cm Long Tube with Silicone lined Straps	500mL	80
UR350050*	LEG BAG - Lever Tap NRV 50cm Long Tube with Silicone lined Straps	500mL	80
UR375010*	LEG BAG - Lever Tap NRV 10cm Short Tube with Silicone lined Straps	750mL	80
UR375030*	LEG BAG - Lever Tap NRV 30cm Long Tube with Silicone lined Straps	750mL	80
UR375050*	LEG BAG - Lever Tap NRV 50cm Long Tube with Silicone lined Straps	750mL	80

\* This product is currently being sold without the tapered shape. All stock will transition to the tapered shape upon depletion of existing stock.

## Leg Bags (T-Tap)

Our Leg Bags are designed to be unobtrusive offering discreet storage under clothing. They have a T-Tap and non-return valve to support a sterile fluid pathway. The soft non-woven backing improves user comfort and they are available with either an adjustable elastic fabric strap or silicone lined strap.

- Discreet to wear
- Soft non-woven backing comfortable next to skin
- Adjustable elastic fabric strap or silicone lined strap
- Multiple tube lengths; 10cm, 30cm and 50cm
- Various volume capacities; 350mL, 500mL and 750mL
- Bonded step connector
- The T-Tap design prevents accidental drainage sometimes experienced with other tap styles



Leg Bags (T-Tap)			STERILE
RE-ORDER	DESCRIPTION	CAPACITY	CARTON QTY
UR135010*	LEG BAG - T-TAP NRV 10cm Short Tube with Silicone Straps	350mL	80
UR135030*	LEG BAG - T-TAP NRV 30cm Long Tube with Silicone Straps	350mL	80
UR135050	LEG BAG - T-TAP NRV 50cm Long Tube with Silicone Straps	350mL	80
UR150010*	LEG BAG - T-TAP NRV 10cm Short Tube with Silicone Straps	500mL	80
UR150030	LEG BAG - T-TAP NRV 30cm Long Tube with Silicone Straps	500mL	80
UR150050*	LEG BAG - T-TAP NRV 50cm Long Tube with Silicone Straps	500mL	80
UR175010	LEG BAG - T-TAP NRV 10cm Short Tube with Silicone Straps	750mL	80
UR175030	LEG BAG - T-TAP NRV 30cm Long Tube with Silicone Straps	750mL	80
UR175050	LEG BAG - T-TAP NRV 50cm Long Tube with Silicone Straps	750mL	80
UR235010*	LEG BAG - T-TAP NRV 10cm Short Tube with Elastic Straps	350mL	100
UR235030*	LEG BAG - T-TAP NRV 30cm Long Tube with Elastic Straps	350mL	100
UR235050*	LEG BAG - T-TAP NRV 50cm Long Tube with Elastic Straps	350mL	100
UR250010	LEG BAG - T-TAP NRV 10cm Short Tube with Elastic Straps	500mL	100
UR250030	LEG BAG - T-TAP NRV 30cm Long Tube with Elastic Straps	500mL	100
UR250050*	LEG BAG - T-TAP NRV 50cm Long Tube with Elastic Straps	500mL	100
UR275010	LEG BAG - T-TAP NRV 10cm Short Tube with Elastic Straps	750mL	100
UR275030*	LEG BAG - T-TAP NRV 30cm Long Tube with Elastic Straps	750mL	100
UR275050*	LEG BAG - T-TAP NRV 50cm Long Tube with Elastic Straps	750mL	100

\* This product is currently being sold without the tapered shape. All stock will transition to the tapered shape upon depletion of existing stock.



## Universal Securement Device

Our Universal Securement Device is designed to optimise the securement of indwelling catheters while ensuring patient comfort. An adhesive strip fastener and soft material tabs provide two levels of securement to ensure the catheter stays in place.

- ⊗ Maximise patient comfort
- ⊗ Two levels of securement to avoid sliding and dislodgement
- ⊗ Low irritant adhesive dressing is soft and comfortable on the skin
- ⊗ Re-adjustable soft material tabs
- ⊗ Simple to apply and remove



## Accessories

Securement of the catheters and drainage bags is a critical aspect of continence care.

We have designed a range of products to simplify securement and ensure bags and catheters perform optimally.

Universal Securement Devices				STERILE
RE-ORDER	DESCRIPTION	DIMENSION	SIZE	CARTON QTY
UR060000	Universal Securement Device	6.5cm x 10.5cm	Standard	50
UR060001	Universal Securement Device	8.5cm x 12.5cm	Large	50



## Urine Bag Hanger

Our easy-to-use reinforced double Urine Bag Hanger features prongs to attach the bag and also includes a clasp to secure the urine bag tubing in place.

- ⊗ Suitable for all major bag types
- ⊗ Reinforced for effective securement
- ⊗ Easy to use prongs to attach the urine bag
- ⊗ Pre-formed securement clasp to secure the urine bag tubing in place



### Urine Bag Hanger

NON-STERILE

RE-ORDER	DESCRIPTION	CARTON QTY
UR063003NS	Urine Bag Hanger	50

## Catheter Tip Syringe

Our Catheter Tip Syringe has a specifically designed funnel end that can be connected to a foley catheter, to irrigate the catheter as clinically indicated.

- ⊗ Ideal for irrigating catheters and other tubes
- ⊗ Clear bold incremental markings
- ⊗ 60mL capacity
- ⊗ Plunger stopper with double sealing ring
- ⊗ Lipid resistant
- ⊗ Free from latex and DEHP



### Catheter Tip Syringe

STERILE

RE-ORDER	DESCRIPTION	VOLUME	CARTON QTY
IV482060	Catheter Tip Syringe	60mL	50





## With you every step of the way.

As an Australian, family owned business, we're committed to supporting the best healthcare outcomes for all Australians, with Continece Care, Venous Access and Airway Management solutions. More than just a manufacturer of devices, we're focused on the end-to-end clinician/patient experience.

From order process, to easy disposal, we do all we can to make our role in the caregiving process seamless.



# 1800 023 321



**m|devices**

27 Llewellyn Avenue  
Villawood NSW 2163  
Australia

T +61 (0)2 8718 2800

F +61 (0)2 8718 2801

[mdevices.com.au](http://mdevices.com.au)

# m|devices®

