



AIRWAY  
MANAGEMENT

Anaesthetic  
& Respiratory  
Products

m|devices®



We believe in effortless performance, so you can focus on caregiving.

One size never fits all. For over a decade we've been designing and manufacturing devices that simplify the roles of clinicians and carers – not just cut to a standard specification but tailored to the unique needs of Australians.

Our focus is on making performance effortless, leaving you free to focus on what matters most – caring for patients.



Leading by Design is a philosophy that reflects our commitment to research and collaboration. A crucial part of our success has been our effort to embed medical professionals in the design process. Their insight helps inform solutions that make the giving of care simple and rewarding.

The problems you face today are the solutions we're working towards for tomorrow.

## Contents

### **Anaesthesia**

#### **Intubation**

- 07 | Laryngeal Masks
- 08 | Laryngeal Masks (Reinforced)
- 09 | Endotracheal Tubes
- 10 | Guedel Airways
- 11 | Intubation Stylets
- 12 | Laryngoscopes Yeescop™
- 13 | Anaesthetic Cushioned Masks
- 14 | Low Profile Anaesthetic Masks

#### **16 | Anaesthetic Circuits**

- 17 | Anaesthetic Circuits

### **Respiratory**

#### **Oxygen Therapy**

- 19 | Oxygen Masks
- 20 | Nasal Oxygen Cannulas
- 21 | Nebuliser Sets
- 22 | Oxygen Tubing
- 23 | Bubble Oxygen Tubing
- 24 | Capnography Masks

#### **Suction Catheters**


- 27 | Suction Catheters Open Tip (Finger Control with Cap)
- 28 | Suction Catheters Open Tip (Y Type Control Vent)
- 29 | Suction Catheters Open Tip (Y Type Control Vent)
- 30 | Suction Catheters Round Tip (Y Type Control Vent)

#### **Suction Tubing**

- 33 | Suction Tubing (Sterile)
- 34 | Suction Tubing (Non-Sterile)

#### **Yankauer Handles**

- 37 | Yankauer Handles



## AIRWAY MANAGEMENT

### Safe, reliable and simple to use

In acute care settings you demand the reliability and performance that comes from solutions carefully tailored to your needs.


Our Airway Management and respiratory support products allow for the simple, safe and reliable achievement and maintenance of airways, the control of fluids and the provision of supplemental oxygen. Our products enable you to maximise patient wellbeing and comfort.



## Laryngeal Mask Airways

A smooth one-piece design allows for easy insertion and guide markings on the handle assist with positioning. The translucent tubing provides clear visibility. The soft-seal cuff reduces pressure and minimises potential trauma. Our masks also feature a 15mm connector and pilot balloon to control the cuff.

- ⊗ One-piece design for easy insertion
- ⊗ Clear guide markers for easy positioning
- ⊗ Soft seal cuff to reduce pressure
- ⊗ Pilot balloon to aid visual and tactile assessment of cuff
- ⊗ 100% medical grade silicone
- ⊗ Translucent tubing for clear visual monitoring
- ⊗ Available in adult, paediatric and neonate sizes
- ⊗ Colour coded for simple size selection



## Intubation

Our full range of intubation devices simplify and reliably support the creation of a viable airway for anaesthesia.

Each product has been carefully considered and designed to aid effective selection, insertion, placement and ongoing management of the airway.


Laryngeal Mask Airways					STERILE
RE-ORDER	DESCRIPTION	SIZE	COLOUR	CARTON QTY	
AN030000	Silicone Disposable Laryngeal Airway Mask	1	Purple	10	
AN030001	Silicone Disposable Laryngeal Airway Mask	1.5	White	10	
AN030002	Silicone Disposable Laryngeal Airway Mask	2	Red	10	
AN030003	Silicone Disposable Laryngeal Airway Mask	2.5	Orange	10	
AN030004	Silicone Disposable Laryngeal Airway Mask	3	Green	10	
AN030005	Silicone Disposable Laryngeal Airway Mask	4	Yellow	10	
AN030006	Silicone Disposable Laryngeal Airway Mask	5	Pink	10	

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Laryngeal Mask Airways (Reinforced)

Our smooth one-piece design allows for easy insertion and the guide markings assist with positioning. Any condensation, blood secretions and regurgitated gastric contents are easily visualised through the translucent tubing. Reinforced tubing reduces kinking and supports more complex intubation.

- ⊗ Reinforced tube for difficult intubations
- ⊗ Available in adult, paediatric and neonate versions
- ⊗ Colour coded for simple size selection
- ⊗ Clear guide markers for easy positioning
- ⊗ Soft seal cuff reduces pressure and minimises potential trauma
- ⊗ A pilot balloon allows for visual and tactile assessment of the inflated / deflated cuff
- ⊗ 100% medical grade silicone
- ⊗ Translucent tubing for clear visual monitoring




Reinforced Laryngeal Mask Airways				STERILE
RE-ORDER	DESCRIPTION	SIZE	COLOUR	CARTON QTY
AN031000	Reinforced Silicone Disposable Laryngeal Airway Mask	1	Purple	10
AN031001	Reinforced Silicone Disposable Laryngeal Airway Mask	1.5	White	10
AN031002	Reinforced Silicone Disposable Laryngeal Airway Mask	2	Red	10
AN031003	Reinforced Silicone Disposable Laryngeal Airway Mask	2.5	Orange	10
AN031004	Reinforced Silicone Disposable Laryngeal Airway Mask	3	Green	10
AN031005	Reinforced Silicone Disposable Laryngeal Airway Mask	4	Yellow	10
AN031006	Reinforced Silicone Disposable Laryngeal Airway Mask	5	Pink	10

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Endotracheal Tubes (Cuffed)

Our Endotracheal Tube has a high volume low pressure cuff. A beveled tip and one-piece design allows for easier insertion through the vocal cords and clear guide markers support an effective positioning. In instances where the tip edge opening makes contact with the tracheal wall and occludes the tube, the Murphy Eye design provides an alternative gas passage.

- ⊗ Large range of sizes
- ⊗ Murphy Eye design provides an alternate gas passage
- ⊗ Radio opaque line for X-ray visualisation
- ⊗ Non-toxic PVC design
- ⊗ High volume cuff ensures positive pressure wall seal
- ⊗ Clear guide markings for easy positioning
- ⊗ Oral/nasal



Endotracheal Tubes (Cuffed)			STERILE
RE-ORDER	DESCRIPTION	SIZE	CARTON QTY
AN021000	Endotracheal Tube - Cuffed	3.0mm	10
AN021001	Endotracheal Tube - Cuffed	3.5mm	10
AN021002	Endotracheal Tube - Cuffed	4.0mm	10
AN021003	Endotracheal Tube - Cuffed	4.5mm	10
AN021004	Endotracheal Tube - Cuffed	5.0mm	10
AN021005	Endotracheal Tube - Cuffed	5.5mm	10
AN021006	Endotracheal Tube - Cuffed	6.0mm	10
AN021007	Endotracheal Tube - Cuffed	6.5mm	10
AN021008	Endotracheal Tube - Cuffed	7.0mm	10
AN021009	Endotracheal Tube - Cuffed	7.5mm	10
AN021010	Endotracheal Tube - Cuffed	8.0mm	10
AN021011	Endotracheal Tube - Cuffed	8.5mm	10
AN021012	Endotracheal Tube - Cuffed	9.0mm	10
AN021013	Endotracheal Tube - Cuffed	9.5mm	10
AN021014	Endotracheal Tube - Cuffed	10mm	10

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Guedel Airways

Our polyethylene, non-sterile oropharyngeal Guedel Airways are colour coded for easy size identification.

- Colour coded to support product selection
- Sized to ISO standards




Guedel Airways				NON-STERILE
RE-ORDER	DESCRIPTION	LENGTH	COLOUR	CARTON QTY
AN090000NS	Guedel Airway - No 000	40mm	Pink	50
AN090001NS	Guedel Airway - No 00	50mm	Blue	50
AN090002NS	Guedel Airway - No 0	60mm	Black	50
AN090003NS	Guedel Airway - No 1	70mm	White	50
AN090004NS	Guedel Airway - No 2	80mm	Green	50
AN090005NS	Guedel Airway - No 3	90mm	Yellow	50
AN090006NS	Guedel Airway - No 4	100mm	Red	50
AN090007NS	Guedel Airway - No 5	110mm	Light Blue	50
AN090008NS	Guedel Airway - No 6	120mm	Orange	50

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Intubation Stylets

Made from malleable aluminium with a medical grade PVC sheath our intubation stylets easily bend to the required shape to assist with tube placement during difficult intubations. A straight tip design assists in tracheal tube exchange.

- Malleable aluminium and PVC sheath design ensures stylets are easily shaped to assist in difficult intubations
- Available in 3 sizes



Intubation Stylets				STERILE
RE-ORDER	DESCRIPTION	SIZE	FOR ETT SIZE	CARTON QTY
AN100000	Intubating Stylet	6Fr	3.0-4.0mm	20
AN100001	Intubating Stylet	10Fr	4.5-7.5mm	20
AN100002	Intubating Stylet	14Fr	8.0-10.0mm	20

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Laryngoscopes - Yeescop™

The Yeescop™ one-piece, single use laryngoscope incorporates handle and blade eliminating handle and blade incompatibility and avoids the Laryngoscope becoming a source of cross contamination.

- ⊗ One-piece design avoids handle and blade incompatibility
- ⊗ Powerful illumination
- ⊗ Single use design avoids cross contamination
- ⊗ Simple function testing



Laryngoscopes - Yeescop™					NON-STERILE
RE-ORDER	DESCRIPTION	SIZE	HANDLE SIZE	BLADE SIZE	CARTON QTY
AN010000	YEESCOPE™ - Disposable Laryngoscope	MAC3	138mm	139mm	10
AN010001	YEESCOPE™ - Disposable Laryngoscope	MAC4	140mm	151mm	10

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Anaesthetic Cushioned Masks

Our Anaesthetic Masks are made from a highly transparent, non-toxic soft PVC ensuring a comfortable facial cushion and an optimal air seal with the patients face. A colour coding system aids product selection and an inflation valve supports efficient and easy inflation and deflation.

- ⊗ Colour coded to aid product selection
- ⊗ Easy grip design
- ⊗ Facial cushion provides a strong seal



Anaesthetic Cushioned Masks				NON-STERILE
RE-ORDER	DESCRIPTION	MASK SIZE	COLOUR	CARTON QTY
AN140001NS	Anaesthesia Mask - Neonate	1	Dark Blue	40
AN140002NS	Anaesthesia Mask - Infant	2	Blue	60
AN140003NS	Anaesthesia Mask - Child	3	Green	50
AN140004NS	Anaesthesia Mask - Small Adult	4	Yellow	40
AN1400045NS	Anaesthesia Mask - Medium Adult	4.5	Orange	40
AN140005NS	Anaesthesia Mask - Adult	5	Pink	40
AN140006NS	Anaesthesia Mask - Large Adult	6	Purple	40

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Low Profile Anaesthetic Masks


Our Low Profile Anaesthetic Masks are made from a highly transparent, non-toxic soft PVC ensuring a comfortable facial fit and an optimal air seal with the patients face. The colour coding system aids product selection and the transparent shell allows visibility.

- ⊗ Colour coded to aid product selection
- ⊗ Easy grip design
- ⊗ Facial cuff provides a strong seal
- ⊗ Transparent shell



## Evolving with your needs.

To deliver on our commitment to make life easier for clinicians, we're constantly developing our offering, taking advantage of the latest technologies and industry best practice.



### Low Profile Anaesthetic Masks

NON-STERILE

RE-ORDER	DESCRIPTION	COLOUR	CARTON QTY
AN145001NS	Anaesthesia Mask - Paediatric	Pink	30
AN145002NS	Anaesthesia Mask - Child	Green	30
AN145003NS	Anaesthesia Mask - Small Adult	Purple	35
AN145045NS	Anaesthesia Mask - Large Adult	Grey	35

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information





## Anaesthetic Circuits

Our circuits are collapsible and are available in two configurations, either corrugated with a rotatable Y-connector or as Bain/Coaxial circuits with a T-connector and a 2L bag.

- ⊗ Two configurations
- ⊗ Available in three lengths
- ⊗ Adult and paediatric versions
- ⊗ Bag limb sold separately



## Anaesthetic Circuits

Our range of circuits have been designed to optimise delivery systems for anaesthesia.

Anaesthetic Circuits are available in a range of different configurations to meet your specific clinical practice needs.

Anaesthetic Circuits			NON-STERILE	
RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY	
AN150003NS	Bain Circuit with T Connector to 60cm Collapsible Circuit	1.8m	20	
AN150004NS	Bain Circuit with T Connector to 60cm Collapsible Circuit	2.4m	20	
AN150005NS	Bain Circuit with T Connector to 60cm Collapsible Circuit	3m	20	
AN150006NS	Corrugated Circuit with a Rotatable Y-Connector	1.8m	20	
AN150007NS	Corrugated Circuit with a Rotatable Y-Connector	2.4m	10	
AN150008NS	Corrugated Circuit with a Rotatable Y-Connector	3m	10	
AN150012NS	Paediatric Circuit with Y-Connector	—	20	
AN151005NS	Bag Limb	3m	10	

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information



## Oxygen Masks

Our Oxygen Masks are made from non-toxic medical grade PVC and are available with and without tubing.

- ⊗ Adjustable strap
- ⊗ Adjustable nose clip
- ⊗ Kink resistant star lumen tubing
- ⊗ Paediatric and adult versions
- ⊗ Up to five-year shelf life
- ⊗ Available in no rebreathing versions with a Reservoir Bag



## Oxygen Therapy

Our Oxygen Therapy solutions are designed to deliver supplemental oxygen across a wide range of clinical settings.

Designed to aid simple and safe administration of oxygen therapy to treat hypoxemia and provide breathing myocardial support.

Oxygen Masks			NON-STERILE
RE-ORDER	DESCRIPTION	SIZE	CARTON QTY
AN06000NS	Oxygen Mask - Standard Shape without Tubing	Paediatric	50
AN060001NS	Oxygen Mask - Standard Shape 2.1m Tubing	Paediatric	50
AN060002NS	Oxygen Mask - Elongated Shape without Tubing	Paediatric	50
AN060003NS	Oxygen Mask - Elongated Shape 2.1m Tubing	Paediatric	50
AN061000NS	Oxygen Mask - Standard Shape without Tubing	Adolescent	50
AN061001NS	Oxygen Mask - Standard Shape 2.1m Tubing	Adolescent	50
AN061002NS	Oxygen Mask - Elongated Shape without Tubing	Adult	50
AN061003NS	Oxygen Mask - Elongated Shape 2.1m Tubing	Adult	50
AN062001NS	Oxygen Mask with Reservoir Bag Elongated Shape 2.1m Tube	Paediatric	50
AN062003NS	Oxygen Mask with Reservoir Bag Elongated Shape 2.1m Tube	Adult	50
AN066001NS	Oxygen Mask - Standard Shape 3m Tubing	Paediatric	50
AN066003NS	Oxygen Mask - Elongated Shape 3m Tubing	Paediatric	50
AN066004NS	Oxygen Mask - Standard Shape 3m Tubing	Adult	50
AN066005NS	Oxygen Mask - Elongated Shape 3m Tubing	Adolescent	50


Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Nasal Oxygen Cannulas

Our Nasal Oxygen Cannulas are made from non-toxic medical grade PVC and are designed to simplify oxygen administration and are soft and flexible to maximise patient comfort.

The oxygen tubing increases visibility and can assist in reducing connection errors.

- ⊗ Kink resistant star lumen
- ⊗ Simple administration
- ⊗ Soft flexible design to maximise patient comfort



Nasal Oxygen Cannulas					NON-STERILE
RE-ORDER	DESCRIPTION	SIZE	TUBE LENGTH	COLOUR	CARTON QTY
AN070000NS	Nasal Oxygen Cannula	Infant	2.1m	Clear	50
AN070001NS	Nasal Oxygen Cannula	Child	2.1m	Clear	50
AN070002NS	Nasal Oxygen Cannula	Adult	2.1m	Clear	50
AN070003NS	Nasal Oxygen Cannula	Infant	3m	Green	50
AN070004NS	Nasal Oxygen Cannula	Child	3m	Green	50
AN070005NS	Nasal Oxygen Cannula	Adult	3m	Green	50

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Nebuliser Sets

All kits are made from non-toxic medical grade PVC and are available with paediatric or adult masks or alternatively with a mouthpiece. Each kit includes a 10mL jar, 2.1m or 3m tubing.

- ⊗ Elongated adult sized masks for better patient comfort
- ⊗ Standard connections
- ⊗ Adjustable strap
- ⊗ Kink resistant star lumen oxygen tubing
- ⊗ Available with mask or mouthpiece
- ⊗ Paediatric and adult versions
- ⊗ Clear volume marking on nebuliser jar



Nebuliser Sets				NON-STERILE	
RE-ORDER	DESCRIPTION	SIZE	TUBE LENGTH	CARTON QTY	
AN080000NS	Nebuliser Kit - Standard Shape 10mL Nebuliser Jar	Paediatric	2.1m	50	
AN080001NS	Nebuliser Kit - Elongated Shape 10mL Nebuliser Jar	Paediatric	2.1m	50	
AN081000NS	Nebuliser Kit - Standard Shape 10mL Nebuliser Jar	Adolescent	2.1m	50	
AN081001NS	Nebuliser Kit - Elongated Shape 10mL Nebuliser Jar	Adult	2.1m	50	
AN082000NS	Nebuliser Kit - Mouth Piece 10mL Nebuliser Jar	-	2.1m	50	
AN083000NS	Nebuliser Kit - Standard Shape 10mL Nebuliser Jar	Paediatric	3m	50	
AN083001NS	Nebuliser Kit - Elongated Shape 10mL Nebuliser Jar	Paediatric	3m	50	
AN084000NS	Nebuliser Kit - Standard Shape 10mL Nebuliser Jar	Adolescent	3m	50	
AN084001NS	Nebuliser Kit - Elongated Shape 10mL Nebuliser Jar	Adult	3m	50	
AN085000NS	Nebuliser Kit - Mouth Piece 10mL Nebuliser Jar	-	3m	50	


Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Oxygen Tubing

Our Oxygen Tubing features a channelled inner lumen to resist kinks and occlusions and is made from non-toxic medical grade PVC.

It is available in four different lengths: 2.1m, 3m, 4m and 30m.

- Soft kink resistant tubing
- Various lengths with optional stepped connector for the 3m length



Oxygen Tubing				NON-STERILE
RE-ORDER	DESCRIPTION	LENGTH	COLOUR	CARTON QTY
AN13000NS	Oxygen Tubing	2.1m	Clear	50
AN130001NS	Oxygen Tubing	3m	Clear	50
AN130004NS	Oxygen Tubing	4m	Clear	50
AN130030NS	Oxygen Tubing	30m	Clear	4
AN135001NS	Oxygen Tubing - Stepped Connector	3m	Clear	50

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Bubble Oxygen Tubing

Our Bubble Oxygen Tubing has incremental sections, varying in diameter along its length. This allows it to be cut to a length which best suits connection to the oxygen source.

- Kink resistant tubing
- Length 30m
- Tubing diameter varies from 3.5mm - 7.5mm



Bubble Oxygen Tubing			NON-STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
AN131030NS	Bubble Oxygen Tubing	30m	4

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Capnography Masks

m|devices Capnography Masks come in a variety of configurations and sizes.

The Adult Elongated Masks include a 10cm extension CO2 sample line outlet and 3m oxygen tubing. There is also the option of a 3m CO2 sample line and /or filter.

In addition to the adult range, a Paediatric Capnography Mask is available. It includes a 10cm extension CO2 sample line outlet, filter and 3m oxygen tubing.

- ⊗ Elongated shape
- ⊗ Adjustable strap
- ⊗ Kink resistant star lumen tubing
- ⊗ Paediatric and adult versions
- ⊗ Up to five-year shelf life
- ⊗ 'All-in-one' product to reduce inventory
- ⊗ Multiple configuration meet clinical needs
- ⊗ Made from non-toxic medical grade PVC



Capnography Masks			NON-STERILE
RE-ORDER	DESCRIPTION	SIZE	CARTON QTY
AN065000NS	Capnography Mask - Filter and 3m Oxygen Tubing	Paediatric	50
AN065001NS	Capnography Mask - Filter and 3m Oxygen Tubing	Adult	50
AN067000NS	Capnography Mask - 3m Oxygen Tubing	Adult	50
AN067001NS	Capnography Mask - 3m CO2 Sample Line, 3m Oxygen Tubing	Adult	50

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information



Designed for medical professionals, by medical professionals.


We know from experience that by observing, listening and collaborating with the right people we can find performance breakthroughs with the potential to transform how you work.



## Suction Catheters Open Tip (Finger Control with Cap)

Designed for easy passage through endotracheal and tracheostomy tubing, our range of Open Tip Catheter with Finger Control and Cap is available in a variety of sizes and offers an open distal end with an elliptical eye for easy and effective suctioning.

- ⊗ Stepped connector
- ⊗ Transparent tubing
- ⊗ Colour coded
- ⊗ Available in a range of sizes
- ⊗ Finger control with cap



Open Tip

## Suction Catheters


Designed to optimise passage through endotracheal and tracheostomy tubes, with minimal irritation to patients, our suction catheters support efficient airway management removing respiratory secretions, preventing plugging and assisting in effective maintenance of the airway.

Suction Catheters Open Tip (Finger Control and Cap)						STERILE
RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON QTY	
AN040000	Suction Catheter, Open Tip, Finger Control and Cap	380mm	5Fr	Grey	50	
AN040001	Suction Catheter, Open Tip, Finger Control and Cap	380mm	6Fr	Light Green	50	
AN040002	Suction Catheter, Open Tip, Finger Control and Cap	560mm	8Fr	Light Blue	50	
AN040003	Suction Catheter, Open Tip, Finger Control and Cap	560mm	10Fr	Black	50	
AN040004	Suction Catheter, Open Tip, Finger Control and Cap	560mm	12Fr	White	50	
AN040005	Suction Catheter, Open Tip, Finger Control and Cap	560mm	14Fr	Green	50	
AN040006	Suction Catheter, Open Tip, Finger Control and Cap	560mm	16Fr	Orange	50	
AN040007	Suction Catheter, Open Tip, Finger Control and Cap	560mm	18Fr	Red	50	

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Suction Catheters Open Tip (Y Type Control Vent)

Designed to assist with a smooth insertion through endotracheal and tracheostomy tubes, for removal of airway secretion/debris that are unable to be spontaneously cleared. The polished open tip range offers soft clear tubing, incremental measurement commencing at 4cm from distal tip. Available in 100mm, 380mm and 560mm lengths to suit your clinical needs.



- ⊗ Polished open tip with bilateral side eyelets
- ⊗ Stepped connector
- ⊗ Soft clear tubing
- ⊗ Colour coded to support product selection
- ⊗ Available in a wide range of diameters
- ⊗ Available in lengths 100mm, 380mm and 560mm
- ⊗ Inner sleeve with perforations to assist with maintaining aseptic technique
- ⊗ Incremental measurement markings




Open Tip

Suction Catheters Open Tip (Y Type Control Vent)						STERILE
RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON QTY	
AN045002	Suction Catheter, Open Tip, Y Type Control Vent	100mm	8Fr	Light Blue	50	
AN045003	Suction Catheter, Open Tip, Y Type Control Vent	100mm	10Fr	Black	50	
AN045004	Suction Catheter, Open Tip, Y Type Control Vent	100mm	12Fr	White	50	
AN042000	Suction Catheter, Open Tip, Y Type Control Vent	380mm	5Fr	Grey	50	
AN042001	Suction Catheter, Open Tip, Y Type Control Vent	380mm	6Fr	Light Green	50	
AN042008	Suction Catheter, Open Tip, Y Type Control Vent	380mm	7Fr	Ivory	50	
AN042002	Suction Catheter, Open Tip, Y Type Control Vent	560mm	8Fr	Light Blue	50	
AN042003	Suction Catheter, Open Tip, Y Type Control Vent	560mm	10Fr	Black	50	
AN042004	Suction Catheter, Open Tip, Y Type Control Vent	560mm	12Fr	White	50	
AN042005	Suction Catheter, Open Tip, Y Type Control Vent	560mm	14Fr	Green	50	
AN042006	Suction Catheter, Open Tip, Y Type Control Vent	560mm	16Fr	Orange	50	
AN042007	Suction Catheter, Open Tip, Y Type Control Vent	560mm	18Fr	Red	50	

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Suction Catheters Open Tip (Y Type Control Vent)

Medium style packaging also available to allow for easy storage in space limited areas.



Open Tip


Suction Catheters Open Tip (Y Type Control Vent)						STERILE
RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON QTY	
AN046000	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	380mm	5Fr	Grey	50	
AN046001	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	380mm	6Fr	Light Green	50	
AN046002	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	8Fr	Light Blue	50	
AN046003	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	10Fr	Black	50	
AN046004	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	12Fr	White	50	
AN046005	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	14Fr	Green	50	
AN046006	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	16Fr	Orange	50	
AN046007	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	18Fr	Red	50	

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Suction Catheters Round Tip (Y Type Control Vent)

Designed to assist with a smooth insertion through endotracheal and tracheostomy tubes. This product comes in three lengths: 100mm, 380mm and 560mm. The round atraumatic tip range offers soft clear tubing, incremental measurement markings commencing at 4cm from distal tip.

- ⊙ Round atraumatic closed tip
- ⊙ Bilateral side eyelets
- ⊙ Available in a wide range of diameters
- ⊙ Available in lengths 100mm, 380mm and 560mm
- ⊙ Inner sleeve with perforations to assist with maintaining aseptic technique
- ⊙ Colour coded for ease of sizing selection
- ⊙ Stepped connector
- ⊙ Transparent
- ⊙ Incremental measurement markings



Round Tip

### Suction Catheters Round Tip (Y Type Control Vent)

STERILE

RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON QTY
AN044002	Suction Catheter, Round Tip, Y Type Control Vent	100mm	8Fr	Light Blue	50
AN044003	Suction Catheter, Round Tip, Y Type Control Vent	100mm	10Fr	Black	50
AN044004	Suction Catheter, Round Tip, Y Type Control Vent	100mm	12Fr	White	50
AN043000	Suction Catheter, Round Tip, Y Type Control Vent	380mm	5Fr	Grey	50
AN043001	Suction Catheter, Round Tip, Y Type Control Vent	380mm	6Fr	Light Green	50
AN043002	Suction Catheter, Round Tip, Y Type Control Ventt	560mm	8Fr	Light Blue	50
AN043003	Suction Catheter, Round Tip, Y Type Control Vent	560mm	10Fr	Black	50
AN043004	Suction Catheter, Round Tip, Y Type Control Vent	560mm	12Fr	White	50
AN043005	Suction Catheter, Round Tip, Y Type Control Vent	560mm	14Fr	Green	50
AN043006	Suction Catheter, Round Tip, Y Type Control Vent	560mm	16Fr	Orange	50
AN043007	Suction Catheter, Round Tip, Y Type Control Vent	560mm	18Fr	Red	50

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

## Where confidence meets care.

Our focus is on making performance effortless, leaving you free to focus on what matters most – caring for patients.








## Suction Tubing (Sterile)

Our sterile Suction Tubing is available in flexible, ribbed and special soft versions and comes in a variety of lengths. Each tube comes with standard hexangular connectors attached.

This product is well suited for surgical use.

- ⊗ Flexible, ribbed and special soft tubing options
- ⊗ Choice of tubing lengths
- ⊗ Hexangular connects
- ⊗ Kink and crush resistant tubing that maintain lumen patency
- ⊗ Securely packaged with double wrap
- ⊗ Single use, sterile and ready to use

Independent testing exceeds minimum requirements for suction tubing degree of collapse as per AS 2120.3 Clauses 5.2.2



## Suction Tubing

Our transparent Suction Tubing is designed to prevent kinking and crushing of tube walls whilst maintaining lumen patency.

Suction Tubing			STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
AN050000	Suction Tubing - Flexible, ID6mm OD9mm	2m	80
AN050001	Suction Tubing - Flexible, ID6mm OD9mm	3m	50
AN050002	Suction Tubing - Flexible, ID6mm OD9mm	4.5m	40
AN050003	Suction Tubing - Flexible, ID6mm OD9mm	6m	35
AN050004	Suction Tubing - Flexible, ID6mm OD9mm	8m	25
AN050009	Suction Tubing - Flexible, ID6mm OD9mm	1m	50
AN051000	Suction Tubing - Ribbed, ID6mm OD9mm	2m	80
AN051001	Suction Tubing - Ribbed, ID6mm OD9mm	3m	50
AN051002	Suction Tubing - Ribbed, ID6mm OD9mm	4.5m	40
AN051003	Suction Tubing - Ribbed, ID6mm OD9mm	6m	35
AN051004	Suction Tubing - Ribbed, ID6mm OD9mm	8m	25
AN054002	Suction Tubing - Ribbed, ID6.8mm OD10mm	2m	80
AN054006	Suction Tubing - Ribbed, ID6.8mm OD10mm	4m	50
AN056001	Suction Tubing - Special Soft, 1 connector	3m	55
AN056006	Suction Tubing - Special Soft, 1 connector	4m	50
AN057001	Suction Tubing - Special Soft, 2 connectors	3m	55
AN057006	Suction Tubing - Special Soft, 2 connectors	4m	50
AN059000	Suction Tubing - Special Soft, with stepped connector	0.48m	120

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information



## Suction Tubing (Non-Sterile)

Our non-sterile Suction Tubing is available in flexible, ribbed and special soft versions and comes in a variety of lengths. Each tube comes with standard hexangular connectors attached.

This product is well-suited for use in non-sterile environments such as wards.

- ⊙ Flexible, ribbed and special soft tubing options
- ⊙ Hexangular connects
- ⊙ Kink and crush resistant tubing that maintain lumen patency
- ⊙ Variable lengths

Independent testing exceeds minimum requirements for suction tubing degree of collapse as per AS 2120.3 Clauses 5.2.2



Suction Tubing			NON STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
AN050000YNS	Suction Tubing - Flexible ID6mm OD9mm	2m	80
AN050001YNS	Suction Tubing - Flexible ID6mm OD9mm	3m	65
AN050002YNS	Suction Tubing - Flexible ID6mm OD9mm	4.5m	45
AN050003YNS	Suction Tubing - Flexible ID6mm OD9mm	6m	35
AN050004YNS	Suction Tubing - Flexible ID6mm OD9mm	8m	30
AN051000YNS	Suction Tubing - Ribbed, ID6mm OD9mm	2m	80
AN051001YNS	Suction Tubing - Ribbed, ID6mm OD9mm	3m	65
AN051002YNS	Suction Tubing - Ribbed, ID6mm OD9mm	4.5m	45
AN051003YNS	Suction Tubing - Ribbed, ID6mm OD9mm	6m	35
AN051004YNS	Suction Tubing - Ribbed, ID6mm OD9mm	8m	30
AN056001YNS	Suction Tubing - Special Soft	3m	65
AN056006YNS	Suction Tubing - Special Soft	4m	50
AN057001YNS	Suction Tubing - Special Soft	3m	50
AN059000NS	Suction Tubing - Special Soft with stepped connector	0.48m	150

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information

We believe  
in effortless  
performance.



By understanding your clinical needs, we've been able to design and develop reliable, high-quality solutions that work effortlessly with your standard of care.



## Yankauer Handles

Our Yankauer Handles are available with straight and crown tips and with and without vents and with an on/off switch. Each comes with a step connector end to ensure a firm connection to any standard suction tubing. They are also available as standard and double wrapped versions.

- ⊗ Smooth and clean edges for atraumatic suctioning
- ⊗ Stepped connectors
- ⊗ Vented and unvented options
- ⊗ Straight and crown tip options
- ⊗ Double wrapped options
- ⊗ Transparent
- ⊗ Paediatric and adult versions



## Yankauer Handles

Designed to allow effective suction without damage to surrounding tissue our Yankauer Handles simply and reliably remove pharyngeal secretions. A transparent handle permits visual monitoring of suction fluids.

Yankauer Suction Handles						
RE-ORDER	DESCRIPTION	OUTER DIAMETER	TIP	STERILE /NS	CARTON QTY	
AN170000	Yankauer Handle - Straight Tip without Vent - Single Wrap	8mm	Straight	Sterile	50	
AN170001	Yankauer Handle - Straight Tip with Vent - Single Wrap	8mm	Straight	Sterile	50	
AN170005	Yankauer Handle - Straight Tip with On/Off Switch - Single Wrap	8mm	Straight	Sterile	50	
AN174000	Yankauer Handle - Straight Tip without Vent - Double Wrap	8mm	Straight	Sterile	50	
AN174001	Yankauer Handle - Straight Tip with Vent - Double Wrap	8mm	Straight	Sterile	50	
AN174005	Yankauer Handle - Straight Tip with On/Off Switch - Double Wrap	8mm	Straight	Sterile	50	
AN170002	Yankauer Handle - Crown Tip without Vent - Single Wrap	8mm	Crown	Sterile	50	
AN170003	Yankauer Handle - Crown Tip with Vent - Single Wrap	8mm	Crown	Sterile	50	
AN170004	Yankauer Handle - Crown Tip with On/Off Switch - Single Wrap	8mm	Crown	Sterile	50	
AN174002	Yankauer Handle - Crown Tip without Vent - Double Wrap	8mm	Crown	Sterile	50	
AN174003	Yankauer Handle - Crown Tip with Vent - Double Wrap	8mm	Crown	Sterile	50	
AN174004	Yankauer Handle - Crown Tip with On/Off Switch - Double Wrap	8mm	Crown	Sterile	50	
AN170007	Yankauer Handle - Paed. Fine Tip with Vent - Single Wrap	4mm	Paed. Fine	Sterile	50	
AN170008	Yankauer Handle - Paed. Fine Tip without Vent - Single Wrap	4mm	Paed. Fine	Sterile	50	
AN170009	Yankauer Handle - Paed. Reg. Tip with Vent - Single Wrap	6mm	Paed. Reg.	Sterile	50	
AN170010	Yankauer Handle - Paed. Reg. Tip without Vent - Single Wrap	6mm	Paed. Reg.	Sterile	50	
AN174007	Yankauer Handle - Paed. Fine Tip with Vent - Double Wrap	4mm	Paed. Fine	Sterile	50	
AN174008	Yankauer Handle - Paed. Fine Tip without Vent - Double Wrap	4mm	Paed. Fine	Sterile	50	
AN174009	Yankauer Handle - Paed. Reg. Tip with Vent - Double Wrap	6mm	Paed. Reg.	Sterile	50	
AN174010	Yankauer Handle - Paed. Reg. Tip without Vent - Double Wrap	6mm	Paed. Reg.	Sterile	50	

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information (Paed. = Paediatric)



## With you every step of the way.

As an Australian, family owned business, we're committed to supporting the best healthcare outcomes for all Australians, with Continece Care, Venous Access and Airway Management solutions. More than just a manufacturer of devices, we're focused on the end-to-end clinician/patient experience.

From order process, to easy disposal, we do all we can to make our role in the caregiving process seamless.

# 1800 023 321



**m|devices**

27 Llewellyn Avenue  
Villawood NSW 2163  
Australia

T +61 (0)2 8718 2800

F +61 (0)2 8718 2801

[mdevices.com.au](http://mdevices.com.au)

# m|devices<sup>®</sup>

m|devices is the trading name of  
Multigate Medical Devices Pty Ltd.  
ABN 69 132 290 058

carbonneutral<sup>®</sup>



A Carbon  
Neutral  
Company

