


Endotracheal Tubes (Cuffed)

Our Endotracheal Tube has a high volume low pressure cuff. A beveled tip and one-piece design allows for easier insertion through the vocal chords and clear guide markers support an effective positioning. In instances where the tip edge opening makes contact with the tracheal wall and occludes the tube, the Murphy Eye design provides an alternate gas passage.

- ⊗ Large range of sizes
- ⊗ Murphy Eye design provides an alternate gas passage
- ⊗ Radio opaque line for X-ray visualisation
- ⊗ Non-toxic PVC design
- ⊗ High volume cuff ensures positive pressure wall seal
- ⊗ Clear guide markings for easy positioning
- ⊗ Oral/nasal



Endotracheal Tubes (Cuffed)			STERILE
RE-ORDER	DESCRIPTION	GAUGE	CARTON QTY
AN021000	Endotracheal Tube - Cuffed	3.0mm	10
AN021001	Endotracheal Tube - Cuffed	3.5mm	10
AN021002	Endotracheal Tube - Cuffed	4.0mm	10

Table continued on next page

Endotracheal Tubes (Cuffed)

Endotracheal Tubes (Cuffed)			STERILE
RE-ORDER	DESCRIPTION	GAUGE	CARTON QTY
AN021003	Endotracheal Tube - Cuffed	4.5mm	10
AN021004	Endotracheal Tube - Cuffed	5.0mm	10
AN021005	Endotracheal Tube - Cuffed	5.5mm	10
AN021006	Endotracheal Tube - Cuffed	6.0mm	10
AN021007	Endotracheal Tube - Cuffed	6.5mm	10
AN021008	Endotracheal Tube - Cuffed	7.0mm	10
AN021009	Endotracheal Tube - Cuffed	7.5mm	10
AN021010	Endotracheal Tube - Cuffed	8.0mm	10
AN021011	Endotracheal Tube - Cuffed	8.5mm	10
AN021012	Endotracheal Tube - Cuffed	9.0mm	10
AN021013	Endotracheal Tube - Cuffed	9.5mm	10
AN021014	Endotracheal Tube - Cuffed	10mm	10