


**Suction Catheters Round Tip (Y Type Control Vent)**

Designed to assist with a smooth insertion through endotracheal and tracheostomy tubes, for the removal of airway secretion/debris that are unable to be spontaneously cleared.

This product comes in three sizes: 100mm, 380mm and 560mm. The round atraumatic tip range offers soft clear tubing, incremental measurement markings commencing at 4cm from distal tip.

- ⊗ Round atraumatic closed tip
- ⊗ Polished round tip with side eyelets
- ⊗ Available in a wide range of diameters
- ⊗ Available in lengths 100mm, 380mm and 560mm
- ⊗ Inner sleeve with perforations to assist with maintaining aseptic technique
- ⊗ Colour coded for ease of sizing selection
- ⊗ Stepped connector
- ⊗ Transparent
- ⊗ Incremental measurement markings



Round Tip

| Suction Catheters Round Tip (Y Type Control Vent) |  |        |       |            | STERILE    |
|---|--|--------|-------|------------|------------|
| RE-ORDER  | DESCRIPTION                                      | LENGTH | GAUGE | COLOUR     | CARTON QTY |
| AN044002  | Suction Catheter, Round Tip, Y Type Control Vent | 100mm  | 8Fr   | Light Blue | 50         |
| AN044003  | Suction Catheter, Round Tip, Y Type Control Vent | 100mm  | 10Fr  | Black      | 50         |
| AN044004  | Suction Catheter, Round Tip, Y Type Control Vent | 100mm  | 12Fr  | White      | 50         |

Table continued on next page

**Suction Catheters Round Tip (Y Type Control Vent)**

| Suction Catheters Round Tip (Y Type Control Vent) |  |        |       |             | STERILE    |
|---|--|--------|-------|-------------|------------|
| RE-ORDER  | DESCRIPTION                                      | LENGTH | GAUGE | COLOUR      | CARTON QTY |
| <b>AN043000*</b>                                  | Suction Catheter, Round Tip, Y Type Control Vent | 380mm  | 5Fr   | Grey        | 50         |
| <b>AN043001*</b>                                  | Suction Catheter, Round Tip, Y Type Control Vent | 380mm  | 6Fr   | Light Green | 50         |
| <b>AN043002</b>                                   | Suction Catheter, Round Tip, Y Type Control Vent | 380mm  | 8Fr   | Light Blue  | 50         |
| <b>AN043003</b>                                   | Suction Catheter, Round Tip, Y Type Control Vent | 560mm  | 10Fr  | Black       | 50         |
| <b>AN043004</b>                                   | Suction Catheter, Round Tip, Y Type Control Vent | 560mm  | 12Fr  | White       | 50         |
| <b>AN043005</b>                                   | Suction Catheter, Round Tip, Y Type Control Vent | 560mm  | 14Fr  | Green       | 50         |
| <b>AN043006</b>                                   | Suction Catheter, Round Tip, Y Type Control Vent | 560mm  | 16Fr  | Orange      | 50         |
| <b>AN043007</b>                                   | Suction Catheter, Round Tip, Y Type Control Vent | 560mm  | 18Fr  | Red         | 50         |

\* 5-6Fr Suction Catheters do not feature the polished eyelets available on larger sizes.