


Oxygen Tubing

Our Oxygen Tubing features a channelled inner lumen to resist kinks and occlusions and is made from non-toxic medical grade PVC.

It is available in four different lengths: 2.1m, 3m, 4m and 30m.

- ⊗ Soft kink resistant tubing
- ⊗ Various lengths with optional stepped connector for the 3m length



Oxygen Tubing				NON-STERILE
RE-ORDER	DESCRIPTION	LENGTH	Colour	CARTON QTY
AN130000NS	Oxygen Tubing	2.1m	Clear	50
AN130001NS	Oxygen Tubing	3m	Clear	50
AN130004NS	Oxygen Tubing	4m	Clear	50
AN130030NS	Oxygen Tubing	30m	Clear	4
AN135001NS	Oxygen Tubing - Stepped Connector	3m	Clear	50