


Minimum Volume Extension Sets

Our Minimum Volume Extension Sets have an inner diameter of 1.5mm and are available in a range of lengths.

- ⊗ Transparent, soft, smooth tubing
- ⊗ 1.5mm inner diameter
- ⊗ Female Luer lock to male Luer lock with rotating collar
- ⊗ Made of non-toxic medical grade PVC



Needle Injection Site



NeutralSite™



Needleless Access Site

Minimum Volume Extension Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	ACCESS	CARTON QTY
IV051002	Minimum Volume Extension Set with Female Luer Lock to Male Luer Lock and Rotating Collar (RC)	600cm	-	50
IV051003	Minimum Volume Extension Set with Injection Site Female Luer Lock to Male Luer Lock & RC	165cm	Standard Injection Site	100
IV051004	Minimum Volume Extension Set with NeutralSite™ Luer Lock to Male Luer Lock and RC	165cm	NeutralSite™	100
IV051006	Minimum Volume Extension Set with NeutralSite™, Non-Return valve, Female Luer Lock to Male Luer Lock and RC	20cm	NeutralSite™	100
IV051010	Minimum Volume Extension Set with Anti-Siphon Valve	150cm	-	100