


## Multi Lumen Standard Bore Extension Sets

Our Multi Lumen Standard Bore Extensions Sets offer both 2 and 3 lumen configurations. They have an inner diameter of 3mm and have slide clamps on each lumen. Standard double and triple Sets have male Luer lock with rotating collar to female luer lock ends. Triple Lumen with 3 x NeutralSite also available, featuring male Luer lock with rotating collar.

- ⊗ Female Luer lock to male Luer lock with rotating collar
- ⊗ Slide clamps
- ⊗ Available with 2 and 3 Lumen
- ⊗ NeutralSite access option
- ⊗ 3mm tubing
- ⊗ Transparent medical grade PVC



NeutralSite™

Multi Lumen Standard Bore Extension Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV052300	Multi Lumen Standard Bore Extension Set with 2 x Lumen	15cm	–	100
IV052400	Multi Lumen Standard Bore Extension Set with 3 x Lumen	15cm	–	100
IV052402	Triple Lumen Standard Bore Extension Set with 3 x NeutralSite™	15cm	NeutralSite™	100