


## Stopcocks (Extension Tubing)

Our Stopcock Extension Sets are made from medical grade, non-toxic, kink resistant 3mm diameter tubing and are available in a range of lengths.

- ⊗ 3-way stopcock incorporating extension tubing
- ⊗ 3mm tubing diameter
- ⊗ Male Luer lock with rotating collar
- ⊗ Two female Luer connectors
- ⊗ Flow indicator arrows



3-Way Stopcocks with Extension Tubing				STERILE
RE-ORDER	DESCRIPTION	LENGTH	COLOUR	CARTON QTY
IV110002	Stopcock Lipid Resistant with 15cm Extension Set	15cm	Blue	50
IV110003	Stopcock Lipid Resistant with 25cm Extension Set	25cm	Blue	50
IV110004	Stopcock Lipid Resistant with 50cm Extension Set	50cm	Blue	50
IV110005	Stopcock Lipid Resistant with 80cm Extension Set	80cm	Blue	50
IV110006	Stopcock Lipid Resistant with 100cm Extension Set	100cm	Blue	50