

Insulin Syringes

Our disposable 1mL and 0.5mL syringes have fixed needles and are available in two needle lengths: 8mm and 13mm and are offered in 27G, 29G and 30G.

The retractable safety syringe version retracts the needle into the barrel after injection to prevent re-use and needle-stick injury.

- ⊗ Orange cap for easy identification
- ⊗ Pre-attached needle
- ⊗ Clear barrel
- ⊗ Bold scale markings

The retractable syringe:

- ⊗ Permanently disabled needle reduces risk of needle stick injury
- ⊗ Eliminates blood splatter
- ⊗ Audible lock
- ⊗ Snap off plunger
- ⊗ Sharps waste reduction



| Insulin Syringes | | | | STERILE |
|------------------|---|--------|-----|------------|
| RE-ORDER | DESCRIPTION | VOLUME | IU | CARTON QTY |
| IV120127 | Insulin Syringe with Standard Fixed Needle 27G x 13mm | 1mL | 100 | 100 |
| IV120129 | Insulin Syringe with Standard Fixed Needle 29G x 13mm | 1mL | 100 | 100 |
| IV120130 | Insulin Syringe with Standard Fixed Needle 30G x 13mm | 1mL | 100 | 100 |
| IV120529 | Insulin Syringe with Standard Fixed Needle 29G x 13mm | 0.5mL | 100 | 100 |
| IV121127 | Insulin Syringe with Standard Fixed Needle 27G x 8mm | 1mL | 100 | 100 |
| IV121129 | Insulin Syringe with Standard Fixed Needle 29G x 8mm | 1mL | 100 | 100 |
| IV121130 | Insulin Syringe with Standard Fixed Needle 30G x 8mm | 1mL | 100 | 100 |

Table continued on next page

Insulin Syringes

| Retractable Insulin Syringes | | | | STERILE |
|------------------------------|---|--------|-----|------------|
| RE-ORDER | DESCRIPTION | VOLUME | IU | CARTON QTY |
| IV158029 | Insulin Syringe with Retractable Safety Needle 29G x 8mm | 0.5mL | 50 | 100 |
| IV158129 | Insulin Syringe with Retractable Safety Needle 29G x 8mm | 1mL | 100 | 100 |
| IV150127 | Insulin Syringe with Retractable Safety Needle 27G x 13mm | 1mL | 100 | 100 |
| IV150129 | Insulin Syringe with Retractable Safety Needle 29G x 13mm | 1mL | 100 | 100 |