

Dual Safety Insulin Pen Needles

Our Dual Safety Insulin Pen Needles provide safety and comfort through their advanced technology.

The user is protected from accidental needlestick by the automatic safety shields, which activate at both ends, after injection.

The broad shield diameter offers greater comfort by reducing pressure on the skin and facilitating penetration.

The mdevices Dual Safety Insulin Pen Needles fit a broad variety of pens including; Sanofi Aventis, Lilly, Novo Nordisk, Owen Mumford and Berlin Chemie.

A choice of needle gauges (30G and 31G) and needle lengths (4mm, 5mm, 6mm and 8mm) are available and are coloured coded for easy identification.

- ⊗ Dual automatic shield lock for safety
- ⊗ Less resistance during needle penetration
- ⊗ Lubricated needle for a smooth injection
- ⊗ Several needle gauges and sizes available
- ⊗ Compatible with a broad variety of insulin pens
- ⊗ Packaged in a sealed sterility barrier



Automatic Safety Shield



Automatic Safety Shield

Dual Safety Insulin Pen Needles

| Dual Safety Insulin Pen Needles | | | | STERILE |
|---------------------------------|--------------------------------|--------------------|--------|------------|
| RE-ORDER | DESCRIPTION | SIZE | COLOUR | CARTON QTY |
| IV512305 | Dual Safety Insulin Pen Needle | 0.3mm (30G) x 5mm | Purple | 100 |
| IV512306 | Dual Safety Insulin Pen Needle | 0.3mm (30G) x 6mm | Orange | 100 |
| IV512308 | Dual Safety Insulin Pen Needle | 0.3mm (30G) x 8mm | Blue | 100 |
| IV512314 | Dual Safety Insulin Pen Needle | 0.25mm (31G) x 4mm | Green | 100 |
| IV512315 | Dual Safety Insulin Pen Needle | 0.25mm (31G) x 5mm | Purple | 100 |
| IV512316 | Dual Safety Insulin Pen Needle | 0.25mm (31G) x 6mm | Orange | 100 |