

Triple Safety Insulin Pen Needles

Our Triple Safety Insulin Pen Needles provide maximum safety and comfort through their advanced technology. The triple safety shield design helps to protect the user from needlestick injury. Automatic safety shields activate at both ends after injection, plus this range features of a pre-covered needle tip which provides an extra level of protection.

The broad shield diameter offers greater comfort by reducing pressure on the skin and facilitating penetration.

The mdevices Triple Safety Insulin Pen Needles fit a broad variety of pens including; Sanofi Aventis, Lilly, Novo Nordisk, Owen Mumford and Berlin Chemie.

A choice of needle gauges (30G and 31G) and needle lengths (4mm, 5mm, 6mm and 8mm) are available and are coloured coded for easy identification.



- ⊗ Dual automatic shield lock for safety
- ⊗ Pre-covered needle tip as additional safety feature
- ⊗ Safety indicator is revealed once safety mechanism is activated
- ⊗ Less resistance during needle penetration
- ⊗ Lubricated needle for a smooth injection
- ⊗ Several needle gauges and sizes available
- ⊗ Packaged in a sealed sterile barrier
- ⊗ Compatible with a broad variety of insulin pens



Pre-covered needle tip



Automatic Safety Shield



Automatic Safety Shield

Triple Safety Insulin Pen Needles

Triple Safety Insulin Pen Needles				STERILE
RE-ORDER	DESCRIPTION	SIZE	COLOUR	CARTON QTY
IV523305	Triple Safety Insulin Pen Needle	0.3mm (30G) x 5mm	Purple	100
IV523306	Triple Safety Insulin Pen Needle	0.3mm (30G) x 6mm	Orange	100
IV523308	Triple Safety Insulin Pen Needle	0.3mm (30G) x 8mm	Blue	100
IV523314	Triple Safety Insulin Pen Needle	0.25mm (31G) x 4mm	Green	100
IV523315	Triple Safety Insulin Pen Needle	0.25mm (31G) x 5mm	Purple	100
IV523316	Triple Safety Insulin Pen Needle	0.25mm (31G) x 6mm	Orange	100