



2-Way Foley Catheters

Approved for both urethral and suprapubic catheterisation, our 2-Way Foley is made from 100% medical grade translucent and flexible silicone. A symmetrical balloon design supports reliable positioning and minimises the potential for ridging.

- ⊗ 100% Silicone
- ⊗ Three lengths available: Paediatric 33cm, Female 23cm and Unisex 45cm
- ⊗ Intended for Urethral and SPC use
- ⊗ Long indwelling time of up to 8 weeks*
- ⊗ Colour coded for easy size identification
- ⊗ Symmetrical balloon ensures reliable positioning
- ⊗ Radiopaque
- ⊗ Standard length codes have printed high contrast incremental markings from 6cm to 30cm



Standard Tip

2-Way Foley Catheters				STERILE	
RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY	
UR010006	2-Way Foley Catheter, Standard Tip, 33cm (Paediatric) with 1.5mL Balloon	6Fr	Pink	10	
UR010000	2-Way Foley Catheter, Standard Tip, 33cm (Paediatric) with 3mL Balloon	8Fr	Light Blue	10	
UR010001	2-Way Foley Catheter, Standard Tip, 33cm (Paediatric) with 3mL Balloon	10Fr	Black	10	
UR010012	2-Way Foley Catheter, Standard Tip, 23cm (Female) with 10mL Balloon	12Fr	White	10	
UR010014	2-Way Foley Catheter, Standard Tip, 23cm (Female) with 10mL Balloon	14Fr	Green	10	
UR010016	2-Way Foley Catheter, Standard Tip, 23cm (Female) with 10mL Balloon	16Fr	Orange	10	

Table continued on next page

2-Way Foley Catheters

2-Way Foley Catheters				STERILE
RE-ORDER	DESCRIPTION	GAUGE	COLOUR	CARTON QTY
UR011000	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	12Fr	White	10
UR011001	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	14Fr	Green	10
UR011002	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	16Fr	Orange	10
UR011003	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	18Fr	Red	10
UR011004	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	20Fr	Yellow	10
UR011005	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	22Fr	Purple	10
UR011006	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	24Fr	Blue	10
UR011007	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 10mL Balloon	26Fr	Pink	10
UR011016	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	16Fr	Orange	10
UR011017	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	18Fr	Red	10
UR011018	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	20Fr	Yellow	10
UR011019	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	22Fr	Purple	10
UR011020	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	24Fr	Blue	10
UR011021	2-Way Foley Catheter, Standard Tip, 45cm (Unisex) with 30mL Balloon	26Fr	Pink	10

* As per laboratory-based testing results and validation against requirements ISO 20696.