


Hydrophilic Coated Nelaton Catheters

Our Hydrophilic catheters are made from non-toxic DEHP free PVC. Their smooth water soluble polymer coating combined with a heat polished distal tip supports a smooth insertion without the need of lubricants.

- ⊗ Available in male, female and paediatric lengths
- ⊗ Heat polished distal end aids insertion
- ⊗ Hydrophilic coating supports lubricant free insertion
- ⊗ Colour coded to support easy size identification
- ⊗ DEHP free



| Hydrophilic Coated Nelaton Catheter | | | | STERILE |
|-------------------------------------|--|-------|------------|------------|
| RE-ORDER | DESCRIPTION | GAUGE | COLOUR | CARTON QTY |
| UR043001 | Hydrophilic Coated Nelaton Catheter, 30cm (Paediatric) | 8Fr | Light Blue | 25 |
| UR043002 | Hydrophilic Coated Nelaton Catheter, 30cm (Paediatric) | 10Fr | Black | 25 |

Table continued on next page

Hydrophilic Coated Nelaton Catheters

| Hydrophilic Coated Nelaton Catheter | | | | STERILE |
|-------------------------------------|--|-------|------------|------------|
| RE-ORDER | DESCRIPTION | GAUGE | COLOUR | CARTON QTY |
| UR044000 | Hydrophilic Coated Nelaton Catheter, 18cm (Female) | 8Fr | Light Blue | 25 |
| UR044001 | Hydrophilic Coated Nelaton Catheter, 18cm (Female) | 10Fr | Black | 25 |
| UR044002 | Hydrophilic Coated Nelaton Catheter, 18cm (Female) | 12Fr | White | 25 |
| UR044003 | Hydrophilic Coated Nelaton Catheter, 18cm (Female) | 14Fr | Green | 25 |
| UR044004 | Hydrophilic Coated Nelaton Catheter, 18cm (Female) | 16Fr | Orange | 25 |
| UR044005 | Hydrophilic Coated Nelaton Catheter, 18cm (Female) | 18Fr | Red | 25 |
| UR045000 | Hydrophilic Coated Nelaton Catheter, 40cm (Male) | 8Fr | Light Blue | 25 |
| UR045001 | Hydrophilic Coated Nelaton Catheter, 40cm (Male) | 10Fr | Black | 25 |
| UR045002 | Hydrophilic Coated Nelaton Catheter, 40cm (Male) | 12Fr | White | 25 |
| UR045003 | Hydrophilic Coated Nelaton Catheter, 40cm (Male) | 14Fr | Green | 25 |
| UR045004 | Hydrophilic Coated Nelaton Catheter, 40cm (Male) | 16Fr | Orange | 25 |
| UR045005 | Hydrophilic Coated Nelaton Catheter, 40cm (Male) | 18Fr | Red | 25 |