

Leg Bags - T-Tap

Our Leg Bags are designed to be unobtrusive offering discreet storage under clothing. They have a T-Tap and non-return valve to support a sterile fluid pathway. The soft non-woven backing improves user comfort and they are available with either an adjustable elastic fabric strap or silicone lined strap.

- ⊗ Discreet to wear
- ⊗ Soft non-woven backing is comfortable next to the skin
- ⊗ Adjustable elastic fabric strap or silicone lined strap
- ⊗ Multiple tube lengths; 10cm, 30cm and 50cm
- ⊗ Various volume capacities; 350mL, 500mL and 750mL
- ⊗ Bonded step connector
- ⊗ The T-Tap design prevents accidental drainage sometimes experienced with other tap styles



T-Tap

Leg Bags			STERILE
RE-ORDER	DESCRIPTION	CAPACITY	CARTON QTY
UR135010*	LEG BAG - T-Tap NRV 10cm Short Tube with Silicone Lined Straps	350mL	80
UR135030*	LEG BAG - T-Tap NRV 30cm Long Tube with Silicone Lined Straps	350mL	80
UR135050	LEG BAG - T-Tap NRV 50cm Long Tube with Silicone Lined Straps	350mL	80

Table continued on next page

Leg Bags - T-Tap

Leg Bags			STERILE
RE-ORDER	DESCRIPTION	CAPACITY	CARTON QTY
UR150010*	LEG BAG - T-Tap NRV 10cm Short Tube with Silicone Lined Straps	500mL	80
UR150030	LEG BAG - T-Tap NRV 30cm Long Tube with Silicone lined Straps	500mL	80
UR150050*	LEG BAG - T-Tap NRV 50cm Long Tube with Silicone lined Straps	500mL	80
UR175010	LEG BAG - T-Tap NRV 10cm Short Tube with Silicone lined Straps	750mL	80
UR175030	LEG BAG - T-Tap NRV 30cm Long Tube with Silicone lined Straps	750mL	80
UR175050	LEG BAG - T-Tap NRV 50cm Long Tube with Silicone Lined Straps	750mL	80
UR235010*	LEG BAG - T-Tap NRV 10cm Short Tube with Elastic Straps	350mL	100
UR235030*	LEG BAG - T-Tap NRV 30cm Long Tube with Elastic Straps	350mL	100
UR235050*	LEG BAG - T-Tap NRV 50cm Long Tube with Elastic Straps	350mL	100
UR250010	LEG BAG - T-Tap NRV 10cm Short Tube with Elastic Straps	500mL	100
UR250030	LEG BAG - T-Tap NRV 30cm Long Tube with Silicone Lined Straps	500mL	100
UR250050*	LEG BAG - T-Tap NRV 50cm Long Tube with Elastic Straps	500mL	100
UR275010	LEG BAG - T-Tap NRV 10cm Short Tube with Elastic Straps	750mL	100
UR275030*	LEG BAG - T-Tap NRV 30cm Long Tube with Elastic Straps	750mL	100
UR275050*	LEG BAG - T-Tap NRV 50cm Long Tube with Elastic Straps	750mL	100

** This product is currently being sold without the tapered shape. All stock will transition to the tapered shape upon depletion of existing stock.*