

# Endotracheal Tube Cuffed

## DESCRIPTION

The m|devices Endotracheal Tube - Cuffed are used to assist patients in establishing and maintaining a patent airway via oral or nasal intubation. The standard 15mm connector allows attachment of a variety of breathing systems and anaesthetic circuits. Our comprehensive size range will suit all ages and sizing is printed on pilot balloon to clearly indicate size whilst insitu. These ET tubes are available with an inflatable balloon (or cuff) to secure the tube in place.

## FEATURES & BENEFITS

- Oral/nasal.
- Transparent, soft and smooth
- Murphy eye design provides an alternate gas passage
- Radio opaque line for X-Ray visualise
- High volume, positive cuff pressure
- One way valve for inflation
- Kink resistant
- Large size range
- Bold gradient markings
- Clear guide markings for easy positioning


## PRODUCT DATA

- Pilot Colour: Blue
- Not made with DEHP
- Not made with Natural Rubber Latex
- Material: PVC
- Sterile EO
- Shelf Life: 3 Years
- ARTG Number: 210994

## STANDARDS

- EN ISO 11135-1:2007 Sterilisation of Health care Products
- EN ISO 10993-1:2009 Biological evaluation of medical devices
- EN ISO 14155:2011 Clinical investigation of medical devices for human subjects - Good clinical practice

## PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Murphy Eye
2	Cuff
3	Valve
4	Pilot Balloon
5	Inflating Tube
6	Connector

# Endotracheal Tube Cuffed

Endotracheal Tube - Cuffed			
RE-ORDER	DESCRIPTION	TUBING LENGTH (m)	CARTON/QTY
AN021000	Endotracheal Tube - Cuffed, 3.0mm	190mm	10
AN021001	Endotracheal Tube - Cuffed, 3.5mm	210mm	10
AN021002	Endotracheal Tube - Cuffed, 4.0mm	230mm	10
AN021003	Endotracheal Tube - Cuffed, 4.5mm	240mm	10
AN021004	Endotracheal Tube - Cuffed, 5.0mm	260mm	10
AN021005	Endotracheal Tube - Cuffed, 5.5mm	270mm	10
AN021006	Endotracheal Tube - Cuffed, 6.0mm	290mm	10
AN021007	Endotracheal Tube - Cuffed, 6.5mm	310mm	10
AN021008	Endotracheal Tube - Cuffed, 7.0mm	320mm	10
AN021009	Endotracheal Tube - Cuffed, 7.5mm	330mm	10
AN021010	Endotracheal Tube - Cuffed, 8.0mm	340mm	10
AN021011	Endotracheal Tube - Cuffed, 8.5mm	345mm	10
AN021012	Endotracheal Tube - Cuffed, 9.0mm	350mm	10
AN021013	Endotracheal Tube - Cuffed, 9.5mm	350mm	10
AN021014	Endotracheal Tube - Cuffed, 10.0mm	350mm	10