

# Anaesthetic Circuit

## Bain/Coaxial Circuit with T Connector

### DESCRIPTION

m|devices Anaesthetic Circuits are available in a range of different configurations to meet clinical practical needs. The collapsible Bain/Coaxial circuits with a T Connector and a large 2L bag with different tubing lengths to meet all your clinical needs.


### FEATURES & BENEFITS

- Three lengths available: 1.8m, 2.4m and 3m
- 2L bag
- T Connector to 60cm collapsible circuit

### PRODUCT DATA

- Not made with DEHP
- Material: PE, EVA
- Shelf Life: 3 Years
- ARTG Number: 269092

### PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	T-Connector
2	Bain Coaxial Circuit 1.8m
3	Limb
4	Reservoir Bag
5	Elbow connector with Port Cap Luer

### Anaesthetic Circuit

RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
<b>AN150003NS</b>	Anaesthetic Circuit Bain/Coaxial Circuit with T Connector	1.8m	20
<b>AN150004NS</b>	Anaesthetic Circuit Bain/Coaxial Circuit with T Connector	2.4m	20
<b>AN150005NS</b>	Anaesthetic Circuit Bain/Coaxial Circuit with T Connector	3m	20
<b>AN151005NS</b>	3m Bag Limb - White	3m	10