

Minimum Volume Extension Set

165cm

DESCRIPTION

The Minimum Volume Extension Set is a microbore tubing with female to male luer lock ends and rotating collar. There is a variety of lengths available. The Minimum Volume Extension Set is a 165cm tubing with female to male luer lock ends with rotating collar. It includes a click clamp and a swabable NeutralSite™.

FEATURES & BENEFITS

- Female to fixed male luer lock with rotating collar
- Click clamp
- NeutralSite™
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter: 1.5mm
- Tubing length: 165cm
- Priming volume: 3mL
- Sterile EO
- The product is MRI and CT suitable
- NeutralSite™ Neutral displacement swabable needleless access site
- Not made with DEHP
- Not made with Natural Rubber Latex
- Material: PVC, PE, PC, POM, ABS
- NeutralSite™ Materials: PP, ABS and silicone
- Shelf Life: 3 Years
- ARTG Number: 354371

ISO STANDARDS

- ISO 8536-4 Infusion equipment for medical use- Part 4:Infusion sets for single use, gravity feed
- ISO 594-1 Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment
 - Part 1: General requirements
- ISO 594-2 Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment
 - Part 2:Lock fittings
- ISO14971 Medical Devices- Application of risk management to medical device
- ISO13485 Medical Devices- Quality management systems – Requirements for regulatory purposes

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Female Luer Lock
3	Tubing
4	Clamp
5	Tubing
6	NeutralSite™
7	Tubing
8	Male Luer Lock
9	Cap

Minimum Volume Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV051004	Minimum Volume Extension Set	165cm	100