

# Minimum Volume Extension Set 20cm

## DESCRIPTION

The Minimum Volume Extension Set is a 20cm tubing with female luer lock to male luer lock and rotating collar. It includes a back-check valve and a swabable NeutralSite™.

## FEATURES & BENEFITS

- Female to male luer lock with rotating collar
- Microbore tubing
- Back-check valve
- NeutralSite™
- Single use only
- Peel open pouch


## PRODUCT DATA

- Tubing diameter: 1.5mm
- Tubing length: 20cm
- Priming volume: 0.35mL
- Sterile EO
- NeutralSite™: Neutral displacement swabable needleless access site
- ARTG Number: 354371
- The product is MRI and CT suitable
- Not Made with Natural Rubber Latex
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Material: ABS, LDPE, HDPE, DEHP Free PVC tubing, PC, Silicone
- NeutralSite™ Materials: PP, ABS and Silicone
- Shelf Life: 3 Years

## ISO STANDARDS

- ISO 8536-4 Infusion equipment for medical use- Part 4:Infusion sets for single use, gravity feed
- ISO 594-1 Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment  
—Part 1: General requirements
- ISO 594-2 Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment  
—Part 2:Lock fittings
- ISO14971 Medical Devices- Application of risk management to medical device
- ISO13485 Medical Devices- Quality management systems – Requirements for regulatory purposes

## PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Female Luer Lock
3	Tubing
4	Back Check Valve
5	NeutralSite™
6	Male Luer Lock
7	Cap

## Minimum Volume Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV051006	Minimum Volume Extension Set	20cm	100