

# Multi Lumen Standard Bore Extension Set

15cm  
2 Lumens

## DESCRIPTION

The m|devices Multi Lumen Standard Bore Extension Set with two Lumens, is 15cm in length. It has a female luer lock and slide clamp on each lumen and a male luer lock connector and rotating collar at the other end.

## FEATURES & BENEFITS

- 2 x Female luer lock connectors
- 2 x slide clamps
- 2 x lumen
- Male luer lock with rotating collar
- Peel open pouch
- Single use only


## PRODUCT DATA

- Tubing diameter: I.D- 3mm O.D- 4mm
- Tubing length: 15cm
- Priming volume (2 lumen): 1.75mL
- Lipid resistant
- Sterile EO
- ART6 Number:354371
- The product is MRI and CT suitable
- Not made with DEHP
- Not made with Natural Rubber Latex
- Material: DEHP free PVC, PE, PC, ABS
- Shelf Life: 3 Years

## ISO STANDARDS

- ISO 13485-2016 Quality management systems. Requirements for regulatory purposes
- ISO 1497-2012 Medical devices. Application of risk management to medical devices
- ISO 11135-2014 Sterilisation of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilisation process for medical devices

## PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Female Luer Lock
3	Slide Clamp
4	Tube
5	Connector
6	Tube
7	Male Luer Lock
8	Cap

## Multi Lumen Standard Bore Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV052300	Multi Lumen Standard Bore Extension Set	15cm	100