




Multi Lumen Standard Bore Extension Set

15cm
3 Lumens

PRODUCT DATA SHEET



DESCRIPTION

The Multi Lumen Standard Bore Extension Set with three Lumens, is 15cm in length. It has female luer lock and slide clamp on each lumen and a male luer lock connector and rotating collar at the other end.


FEATURES AND BENEFITS

- 3 x Female luer lock connectors
- 3 x slide clamps
- 3 x lumen
- Male luer lock with rotating collar
- Peel open pouch
- Single use only

PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 15cm
- Priming volume (3 lumen): 2.12mL
- DEHP free
- Sterile EO
- Latex free
- Material: PVC, PE, PC
- Shelf Life: 3 Years
- ARTG Number: 228230
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV052400	Multi Lumen Standard Bore Extension Set	15cm	100