

Triple Lumen Standard Bore Extension Set

3 Lumens
15cm

DESCRIPTION

The Triple Lumen Standard Bore Extension Set enhances your ability to manage fluid deliveries simultaneously. This Triple Lumen set measures 15cm in length and features three NeutralSite™ for needle-free access and a slide clamp on each lumen for precise flow control. For secure connection to access devices, the extension set features a male luer lock connector complete with a rotating collar.

FEATURES & BENEFITS

- 3 NeutralSite™ needleless access sites
- 3 lumens
- 3 slide clamps
- Male luer lock with rotating collar
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 15cm
- Priming volume: 2.1ml
- NeutralSite™: Neutral displacement swabable needleless access site
- Lipid resistant: Yes
- Not Made with DEHP
- Not Made with Natural Rubber Latex
- Material: PVC, ABS, PE, PP, Silicone
- Sterile E0
- Shelf Life: 3 years
- ARTG: 354371

ISO STANDARDS

- ISO 80369-7: Small-bore connectors for liquids and gases in healthcare applications-Part 7: Connectors for intravascular or hypodermic applications
- ISO14971: Medical Devices- Application of risk management to medical device
- ISO13485: Medical Devices- Quality management systems – Requirements for regulatory purposes

PRODUCT DIAGRAM



ITEM	DESCRIPTION		
1	NeutralSite™		
2	Tubing	5	Tubing
3	Slide Clamp	6	Male Luer Lock
4	Connector	7	Cap

Triple Lumen Standard Bore Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV052402	Triple Lumen Standard Bore Extension Set, 3 lumens, 3 NeutralSite™	15cm	100