

Vented Spike Adaptor

10cm with Needleless Access Site

DESCRIPTION

Designed for safe and efficient fluid transfer, the m|devices Vented Spike Adaptor enables the transfer of fluids from bags, vials or bottles without the need for needles. The needle-free access site minimises the risk of the needlestick injuries, promoting user safety.

FEATURES & BENEFITS

- Vented spike
- Standard bore tubing 10cm in length
- Click clamp
- Needleless access site
- Sterile peel open pouch


PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 10cm
- Negative pressure needle free valve
- Sterile EO
- ARTG Number: 354371
- Not Made with Natural Rubber Latex
- Not Made with DEHP
- Material: PE, ABS, PVC, POM, PC, Silicone
- Shelf Life: 3 Years

ISO STANDARDS

- ISO 8536-4:2020 Infusion equipment for medical use- Part 4: Infusion sets for single use, gravity feed
- ISO 11135:2014 Sterilisation of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilisation process for medical devices
- ISO 14971:2019 Medical devices - Application of risk management to medical devices.

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Spike Chamber
2	Tubing
3	Clamp
4	Needleless Access Site

Vented Spike Adaptor			STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV070004	Vented Spike Adaptor	10cm	120