

Cervical Dilating Balloon Procedure Pack

DESCRIPTION

m|devices Cervical Dilating Balloon Procedure Pack has been specifically designed with clinicians to support the insertion of the m|devices Cervical Dilating Balloon. The sterile pack provides all of the necessary items used for this procedure, from the sterile field for the setup trolley through to the inclusion of the m|devices Cervical Dilating Balloon, assisting the clinician with maintaining infection control principles, eliminating collection times and significantly reducing individual packaging waste.

The m|devices Cervical Dilating Balloon features a round closed distal end which meets infection control principles by maintaining a closed system and clear markings for balloon capacity on the balloon lumen and printed batch have all been considered to assist the clinician with ease of use. High contrast position markings have been printed on the shaft to assist clinicians with the placement of the Cervical Dilating Balloon.

FEATURES & BENEFITS

- Specifically designed "fit for purpose" procedure pack
- Inclusion of the m|devices Cervical Dilating Balloon 100% Silicone, closed distal end with high contrast position markings printed on the shaft
- Significant reduction in individual packaging waste
- Pack allows for an efficient setup time and assists with maintaining infection control principles
- Batch tracking sticker provided
- Clinician's valuable time no longer spent collecting individual items

PRODUCT DATA

Contents:

- 1 x Sterile Field 55gsm 120x120cm: CPE
- 1 x Blue rectangle tray:PS
- 2 x 2.7g Lubricating gel sterile: Glycerol and Sodium hydroxide
- 4 x 20mL Syringe Luer Lock: PP
- 1 x 25cm Sponge Forceps (Ramplly): Stainless steel
- 2 x 120mL Gallipots: Modified PP
- 1 x m|devices Cervical Dilating Balloon 80mL 20Fr 40cm: Silicone and PVC NRV
- 1 x XRD Gauze Swab 10x10cm 8 ply (5): cotton with X-Ray detectable filaments
- 1 x 15cm Clear Plastic Forceps: ABS
- 1 x Securement Device 8.5x 12.5cm: Non-woven fabric, velcro, acrylic adhesive and glassine paper
- 1 x Pouch insert: Paper


- ARTG Number: 324486
- Shelf life: 3 Years
- Not Made with Natural Rubber Latex
- EO Sterilisation

STANDARDS

- ISO 14971:2019 Medical devices — Application of risk management to medical devices
- ISO 11135:2014 Sterilisation of health-care products EO

PRODUCT DIAGRAM

Pack Items					
QUANTITY	DESCRIPTION	REF.	QUANTITY	DESCRIPTION	REF.
1	Sterile Field 120x120cm	A	1	m devices Cervical Dilating Balloon	G
1	Blue Rectangle Tray	B	1	XRD Gauze Swab 10x10cm 8Ply (5)	H
2	Lubricating Gel Sterile 27g	C	1	m devices Securement Device	I
4	Luer Lock Syringe 20mL	D	1	Plastic Forceps Clear	J
2	Gallipots 120mL	E	1	Tracking Label	K
1	SS Sponge Forceps 25cm	F			



Cervical Dilating Balloon Procedure Pack

Cervical Dilating Balloon Procedure Pack

RE-ORDER	DESCRIPTION	CARTON/QTY
OB010002	Cervical Dilating Balloon Procedure Pack	30