

Surgical Marker Pens

DESCRIPTION

m|devices Surgical Marking Pens for marking the patient's skin for surgical site selection for diagnostic and interventional radiology, surgery and treatment of skin diseases. The product comes with a ruler and labels and there are two tip size options to meet all clinical requirements.

FEATURES & BENEFITS

- Available as:
Regular tip (1mm) or
Fine tip (0.5mm)
- Long lasting
- Scrub resistant
- Non toxic, non-allergenic and non irritating
- Odourless


PRODUCT DATA

- Radiation sterilised
- Sterile and non-sterile versions
- Ink: Gentian Violet
- Marker: Polypropylene
- ARTG: 210993
- Not Made with Natural Rubber Latex
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Shelf Life: 3 years

STANDARDS

TBC

PRODUCT DIAGRAM



Skin Marking Pen

RE-ORDER	DESCRIPTION	LINE WIDTH	CARTON/QTY
SU010111	Skin Marking Pen with Regular Tip	1.0mm	20
SU010111NS	Skin Marking Pen with Regular Tip - Non Sterile	1.0mm	400
SU010121	Skin Marking Pen with Fine Tip	0.5mm	25