

2-Way Foley Catheter Silicone

DESCRIPTION

Our 2-Way Foley Catheter is made from 100% medical grade translucent and flexible silicone. Unisex and Female lengths are indicated for both urethral and suprapubic catheterisation and our Paediatric sizing is indicated for urethral catheterisation only. A symmetrical balloon design supports reliable positioning and minimises the potential for ridging.

FEATURES & BENEFITS

- Unisex and Female lengths Indicated for both urethral and suprapubic catheterisation
- Paediatric lengths indicated for Urethral catheterisation
- Suitable for paediatric, male and female patients
- Symmetrical balloon
- Smooth round, tapered tip
- Bilateral eyelets for effective drainage
- Radiopaque
- Colour-coded
- Long indwelling time of 8 weeks
- Three lengths available: Paediatric 33cm, Female 23cm and Unisex 45cm
- Unisex length codes have printed high contrast incremental markings from 6cm to 30cm


PRODUCT DATA

- Not Made with Natural Rubber Latex
- Not made with DEHP
- Materials: Silicone
- Shelf life: 3 years
- ARTG number: 211796
- Sterile EO

ISO STANDARDS

- ISO 20696:2018 Sterile urethral catheters for single use
- ISO 14971:2007 Medical devices-Application of risk management to medical devices
- ISO 11135-1:2007 Health-care Product Sterilisation Confirmation and Sterilisation by medical ETO verification and routine control

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Tip
2	Balloon
3	Funnel
4	Colour Specification
5	Side Eyelet

2-Way Foley Catheter Silicone

2-Way Foley Catheter							STERILE
RE-ORDER	DESCRIPTION	SIZE	LENGTH	BALLOON CAPACITY	GAUGE	COLOUR	CARTON/QTY
UR010006	2-Way Foley Catheter	Paediatric	33cm	1.5mL	6Fr	Pink	10
UR010000	2-Way Foley Catheter	Paediatric	33cm	3mL	8Fr	Light Blue	10
UR010001	2-Way Foley Catheter	Paediatric	33cm	3mL	10Fr	Black	10
UR010012	2-Way Foley Catheter	Female	23cm	10mL	12Fr	White	10
UR010014	2-Way Foley Catheter	Female	23cm	10mL	14Fr	Green	10
UR010016	2-Way Foley Catheter	Female	23cm	10mL	16Fr	Orange	10
UR011000	2-Way Foley Catheter	Unisex	45cm	10mL	12Fr	White	10
UR011001	2-Way Foley Catheter	Unisex	45cm	10mL	14Fr	Green	10
UR011002	2-Way Foley Catheter	Unisex	45cm	10mL	16Fr	Orange	10
UR011003	2-Way Foley Catheter	Unisex	45cm	10mL	18Fr	Red	10
UR011004	2-Way Foley Catheter	Unisex	45cm	10mL	20Fr	Yellow	10
UR011005	2-Way Foley Catheter	Unisex	45cm	10mL	22Fr	Purple	10
UR011006	2-Way Foley Catheter	Unisex	45cm	10mL	24Fr	Blue	10
UR011007	2-Way Foley Catheter	Unisex	45cm	10mL	26Fr	Pink	10
UR011016	2-Way Foley Catheter	Unisex	45cm	30mL	16Fr	Orange	10
UR011017	2-Way Foley Catheter	Unisex	45cm	30mL	18Fr	Red	10
UR011018	2-Way Foley Catheter	Unisex	45cm	30mL	20Fr	Yellow	10
UR011019	2-Way Foley Catheter	Unisex	45cm	30mL	22Fr	Purple	10
UR011020	2-Way Foley Catheter	Unisex	45cm	30mL	24Fr	Blue	10
UR011021	2-Way Foley Catheter	Unisex	45cm	30mL	26Fr	Pink	10