

Urine Bag 4000mL

T-Tap Non-Return Valve
150cm Tubing
Integrated Hanger with Rope

PRODUCT DATA SHEET



DESCRIPTION

m|devices Urine Bag has a large 4000mL volume capacity, non-return valve, 150cm tubing length with bonded stepped connector and integrated hanger with rope.

The urine bag is used to collect and contain urine via a urinary catheter or connected leg bag. This sterile bag is ideal for overnight drainage due to its large volume capacity.

The T-Tap outlet allows for ease of use to empty the drainage bag and the clear incremental measurements allow for visual urine output.


FEATURES AND BENEFITS

- 4000mL Large volume capacity
- Long 150cm large Bore Tubing with bonded stepped connector
- Integrated Hanger with Rope for easy attachment to desired location
- Non-return Valve
- Easy to use T-Tap outlet for emptying the bag
- Clear printed incremental markings on bag
- Ideal for overnight drainage

PRODUCT DATA

- Volume: 4000mL
- Tubing Length: 150cm
- Tubing ID/OD: 9 x 1500mm (inlet tubing)
- Material: PVC, PP, PE, Polyester rope, Nylon air filter
- Latex free
- Sterile EO
- DEHP: Yes
- ARTG Number: 267875
- Shelf life: 5 Years

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	LENGTH	BAG CAPACITY	CARTON/QTY
UR064000	Urine Bag T-Tap NRV Integrated Hanger with Rope	150cm	4000mL	40pcs